

PHILIPS

EnergyLight



HF3309/01



Feel as energetic

as on a sunny day

Special light device with intensity dimmer that gives you the daylight—and energy—you need, whatever the season, in a way that easily fits into your daily routine.

Feel energetic in 30 minutes a day

- Original Philips EnergyLight lamps

Digital feedback on increased energy

- Digital timer and display

Comfortable light

- Light intensity dimmer

Easy positioning for personal preference


- Adjustable angle

Safe to use

- UV-free

Highlights

Original Philips EnergyLight lamps



The EnergyLight features original Philips EnergyLight lamps producing light with all the intensity and quality of natural daylight. In fact, it can give you up to 10,000 lux, but without any unpleasant dazzle or glare. If you use the EnergyLight for just 30 minutes a day (there's no need to stare into the lamp all the time) at a light intensity of 10,000 lux, you'll notice an improvement in your energy level and mood after just 5 days' use.

Digital timer and display



The EnergyLight HF 3309 is equipped with a digital user interface with built-in user protocol for easy use, which also gives direct feedback on the amount of energy gained.

Light intensity dimmer



Fight winter blues and get the energy you have in summer while using the appliance at a convenient level of light you need for easy-

reading etc. The EnergyLight features a light intensity dimmer so it can easily be adjusted between zero to about 10,000 lux (at 15-20 cm) matching the intensity of your personal preference.

Adjustable angle



The appliance can be positioned at various angles for comfortable use.

UV-free



The EnergyLight uses the full natural daylight spectrum. It meets European safety standards and features a UV-filtering screen. The EnergyLight is therefore safe to use.

Supported by the Medical Device Directive (MDD)



This device has proven health benefits, as certified by the Medical Device Directive (MDD).

Specifications

Technical specifications

Power: 72 W

Voltage: 220/230

Frequency: 50/60 Hz

Type of lamps: 2 PL-L 36W Philips EnergyLight

Lifetime of lamps: 10.000

Type of ballast: HF-R 232PLT (High Frequency)

Ballast operating frequency: 42.000 kHz

Insulation: Class II (double isolation)

Cord length: 3 m

Weight and dimensions

Product dimensions: 33x57x18 (WxHxD) cm

Product weight: 3.5 kg

F-box dimensions (= A-box): 36x62.2x21.5 (WxHxD) cm

F-box weight (= A-box): 4.4 kg

Qnt. on Euro pallet: 30 (UK 39) pcs

Logistic data

CTV code: 884330901000

Country of origin: Germany

Safety

UV-free: No UV comes through the filtering screen

No flickering: due to high frequency ballast avoiding headaches

Cenelec certified

Easy to use

On product user protocol

Digital user interface: for feedback on energy and time

Energybar: indicating how much energy is gained

Digital timer

Pause function

On/off button

Feel energetic in 30 min a day

Original Philips EnergyLight: lamps: 10.000 lux at 20 cm, 2.500 Lux at 60 cm

Comfortable light

Light intensity dimmer

Easy positioning

Adjustable angle: 20 degrees

Compact size: 18x33 (DxW) cm

Satisfaction indicator

94% satisfied users

Best buy: award by STIWA (2003) as successor of HF 3305

Medical appliance

Proven health benefits

