

PHILIPS

sonicare

Sonic electric
toothbrush with app

DiamondClean
Smart

Personalised coaching

Smart brush head sensors

Smart brush head recognition

5 modes, 3 intensities



HX9924/14



Complete system for a healthier mouth

Our best toothbrush with smart sensor technology

Get the ultimate clean with DiamondClean Smart electric toothbrush. Eliminate guesswork with real-time feedback and personalised coaching for superior teeth, gums, whitening and fresh breath. It's more than just a toothbrush - it's a complete oral care solution.

Proven to improve oral health

- Removes up to 10 x more plaque* for a deep clean
- Up to 7 x healthier gums* in just 2 weeks
- Remove up to 100% more stains** in just 3 days
- Enjoy long-lasting fresh breath

Personalised in-app coaching to perfect brushing

- Guides your brushing and highlights missed areas
- Start taking better care of your problem areas
- Alerts you if you're pressing too hard
- Guides you to reduce your scrubbing
- Takes you back to any missed spots for 100% coverage

Optimise your brushing

- Brush heads automatically select the optimal settings
- Make sure you're getting the most from your brush head
- Choose from 5 modes and 3 intensity settings



asimpleswitch.com


Highlights

Say goodbye to plaque



Click on the C3 Premium Plaque Control brush head to experience our deepest clean. With soft flexible sides, the bristles contour to the shape of each tooth to provide 4 x more surface contact*** and remove up to 10 x more plaque from hard-to-reach areas*.

Care for your gums



Click on the G3 Premium Gum Care brush head to improve your gum health. Its smaller size and targeted gumline bristles deliver gentle yet effective cleaning along the gum line, where gum disease starts. It's clinically proven to deliver up to 100% less gum inflammation* and up to 7 x healthier gums in just two weeks.*

Whiter teeth, fast




Click on the W3 Premium White brush head to remove surface stains and reveal your whitest smile. With its densely packed central stain removal bristles, it's clinically proven to remove up to 100% more stains in just three days.**

TongueCare+ tongue brush




Click on the TongueCare+ tongue brush to gently remove odour-causing bacteria from the pores of your tongue. Its 240 specially designed MicroBristles get between all your tongue's ridges and grooves to remove bacteria and debris that lead to bad breath. Team with our antibacterial BreathRx tongue spray for superior cleaning and super-fresh breath.

Never miss a spot




The Philips Sonicare app lets you know when you've achieved a thorough clean and coaches you to be a better, more mindful brusher.

Focus areas and goal-setting



Any problem areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map and remind you to pay extra attention to these areas.

Too much pressure?



You might not notice when you brush too hard, but your DiamondClean Smart will. If you use too much pressure, the light ring on the end of your handle will flash. This is a gentle reminder to ease off the pressure and let your brush head do the work. 7 out of 10 people found that this feature helped them to become a better brusher.

Easy does it



With an electric toothbrush, you let the brush do the work. We've built a scrubbing sensor into your handle as a gentle reminder to reduce scrubbing. This way, you can improve your technique and get a gentler, more effective clean.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Modes

Clean: For exceptional everyday cleaning

White+: To remove surface stains

Deep Clean+: For an invigorating deep clean

Gum Health: Special attention to molars

TongueCare+: For long-lasting clean breath

3 intensities: High, Medium, Low

Items included

Handle: 1 DiamondClean Smart

Brush heads: 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care, 1 W3 Premium White, 1 TongueCare+ tongue brush

Glass charger: 1

Travel case: USB charging travel case

Antibacterial tongue spray: 2 oz

BrushSync Mode Pairing

C3 Premium Plaque Control: Pairs with Clean mode

G3 Premium Gum Care brush head: Pairs with Gum Health mode

W3 Premium White brush head: Pairs with White+ mode

TongueCare+ tongue brush: Pairs with TongueCare mode

Connectivity

Bluetooth® wireless technology: Connected brushing app

Compatibility

Android compatibility: Android phones, Bluetooth 4.0-enabled tablets

iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher

Design and finishing

Colour: Black

Cleaning performance

Speed: Up to 62,000 brush movements/min

Performance: Removes up to 10 x more plaque*

Health benefits: Up to 7 x healthier gums****

Whitening benefits: Up to 100% less stains*****

Pressure feedback: Ring lights up in purple, Vibrates handle to alert user

TouchUp: Ensures 100% coverage

Timer: BrushPacer and Smart-timer

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

Location sensor: Tracks and improves coverage

Scrubbing sensor: Guides to reduce scrubbing

3D mouth map: Highlights trouble areas

TouchUp: To address missed spots

Position sensor: Guides to areas you missed

Ease of use

Handle compatibility: Easy click-on brush heads

Brush head recognition: Syncs with optimal mode

Replacement reminder: To always ensure best results, reminder icon lights up

Battery indicator: Illuminated icon indicates battery life

Handle: Slim ergonomic design

Display: Illuminated display

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): 14 days*****

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

Ingredients

Antibacterial tongue spray: 0.09% CPC, Zytex



* Removes up to 7x more plaque vs. a manual toothbrush

* * as compared to a manual toothbrush with a leading whitening toothpaste

* * compared to DiamondClean

* * * in Gum Health mode vs a manual toothbrush at 2 weeks

* * * * after 3 days of use in White+ mode vs a manual toothbrush

* * * * based on two periods of two-minute brushings per day on standard mode