

PHILIPS

sonicare

Sonic electric
toothbrush with app

DiamondClean
Smart

Personalized coaching

Smart brush head sensors

Smart brush head recognition

5 modes, 3 intensities



HX9924/62



Complete system for a healthier mouth

Our best toothbrush with smart sensor technology

Get the ultimate clean with DiamondClean Smart. Eliminate guesswork with real time feedback and personalized coaching for superior teeth, gum, whitening & fresh breath. It's more than just a toothbrush - it's a complete oral care solution.

Proven to improve oral health

- All-in-One brush head for complete care
- Up to 7x healthier gums**** in just 2 weeks
- Enjoy long-lasting fresh breath

Personalized in-app coaching to perfect brushing

- Start taking better care of your trouble areas
- Guides you to brush and highlights missed areas
- Alerts you if you're pressing too hard
- Guides you to reduce your scrubbing
- Takes you back to any missed spots for 100% coverage

Optimize your brushing

- Brush heads automatically select the optimal settings
- Make sure you're getting the most from your brush head
- Choose from 5 modes and 3 intensity settings
- Premium charging glass and travel case


Highlights

Premium All-in-One brush head



Clean without compromise with the A3 Premium All-in-One brush head. Our all-time best has angled bristles to help remove up to 20x more plaque even in hard-to-reach areas*. Triangular tips work to remove up to 100% more stains in less than two days***. And longer bristles clean deep for gums that are up to 15x healthier in just two weeks**. All that without ever switching brush heads. Did you know: brush heads become less effective after 3 months of use? Our BrushSync™ feature can remind you when it's almost time for a replacement.

Care for your gums




Click on the G3 Premium Gum Care brush head to improve your gum health. Its smaller size and targeted gum line bristles deliver gentle yet effective cleaning along the gum line, where gum disease starts. It's clinically proven to deliver up to 100% less gum inflammation**** and up to 7x healthier gums in just two weeks.****

TongueCare+ tongue brush




Click on the TongueCare+ tongue brush to gently remove odor causing bacteria from the pores of your tongue. Its 240 specially designed MicroBristles get between all your tongue's ridges and grooves to remove bacteria and debris that lead to bad breath. Team with our antibacterial BreathRx tongue spray for superior cleaning and super fresh breath.

Focus areas and goal-setting




Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map, and remind you to pay extra attention to these areas.

Never miss a spot



The Philips Sonicare app lets you know when you've achieved a thorough clean, and coaches you to be a better, more mindful brusher.

Too much pressure?




You might not notice when you brush too hard, but your DiamondClean Smart will. If you use too much pressure, the light ring on the end of your handle will flash. This is a gentle reminder to ease off the pressure, and let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

Easy does it



With an electric toothbrush, you let the brush do the work. We've built a scrubbing sensor into your handle as a gentle reminder to reduce scrubbing. This way, you can improve your technique and get a gentler, more effective clean.

TouchUp feature



If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.

Specifications

Modes

Clean: For exceptional everyday clean

White+: To remove surface stains

Deep Clean+: For an invigorating deep clean

Gum Health: Special attention to molars

TongueCare+: For long lasting clean breath

3 intensities: Low, Medium, High

Items included

Handle: 1 DiamondClean Smart

Brush heads: 2 A3 Premium All-in-One, 1 G3 Premium Gum Care, 1 TongueCare+ tongue brush

Travel case: Premium travel case

Glass charger: 1

BrushSync Mode Pairing

C3 Premium Plaque Control: Pairs with Clean mode

G3 Premium Gum Care brush head: Pairs with Gum Health mode

W3 Premium White brush head: Pairs with White+ mode

TongueCare+ tongue brush: Pairs with TongueCare mode

Connectivity

Bluetooth® wireless technology: Connected brushing app

Compatibility

Android compatibility: Android phones, Bluetooth 4.0 enabled tablets

iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher

Design and finishing

Color: Rose gold

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Removes up to 20x more plaque*

Health benefits: Up to 15x healthier gums**

Whitening benefits: Up to 100% less stains***

Pressure feedback: Ring lights up in purple, Vibrates handle to alert user

TouchUp: Ensures 100% coverage

Timer: BrushPacer and SmarTimer

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

Location sensor: Tracks & improves coverage

Scrubbing sensor: Guides to reduce scrubbing

3D mouth map: Highlights trouble areas

TouchUp: To address missed spots

Ease of use

Handle compatibility: Easy click-on brush heads

Brush head recognition: Syncs with optimal mode

Replacement reminder: To always ensure best results, reminder icon lights up

Battery indicator: Illuminated icon indicates battery life

Handle: Slim ergonomic design

Display: Illuminated display

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): 14 days****

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty



* vs. a manual toothbrush

** vs. a manual toothbrush in 6 weeks

*** in lab test vs. a manual toothbrush

**** in Gum Health mode vs a manual toothbrush at 2 weeks

***** based on two periods of two-minute brushings per day on standard mode