

# PHILIPS

## sonicare

Sonic electric  
toothbrush

DiamondClean

2 brush heads

Glass charger

with Polish mode



HX9352/49



## Whiter, healthier teeth for life

### Best Philips Sonicare whitening sonic toothbrush

Philips Sonicare's best whitening in our most elegant sonic electric toothbrush. Switch to Philips Sonicare.

#### Provides a superior clean

- Up to 10x more plaque removal than a manual toothbrush

#### Proven to improve oral health

- Whitens teeth in just one week
- Improves gum health in just 2 weeks\*
- Choose from 5 modes including Polish and Sensitive

#### Innovative technology

- Philips Sonicare's advanced sonic technology

#### Optimize your brushing

- Safe & gentle on sensitive areas, orthodontics & dental work

#### Designed around you


- Easy start program builds your Philips Sonicare routine
- Timers encourage thorough 2 minutes of brushing
- USB charging travel case and charging glass



[asimpleswitch.com](http://asimpleswitch.com)


# Highlights

## Say goodbye to plaque




And thanks to its flexible design, it removes up to 10x more plaque from hard-to-reach areas\* for a truly deep clean along the gum line and between teeth.

## Look after those gums



With optimal cleaning from your W DiamondClean, gums will become healthier in 2 weeks\*. Removing up to 7x more plaque along the gum line than a manual toothbrush\*, you'll get your healthiest smile.

## Naturally whiter teeth



Click on the W2 Optimal White brush head to remove surface stains and reveal a whiter smile. With its densely packed central stainremoval bristles, it's clinically proven to whiten teeth in just one week.

## 5 modes to clean your way




With DiamondClean, you get a refreshing clean every day. Our 5 modes cover all your brushing needs; Clean mode – for exceptional daily cleaning, Gum Care – to gently massage gums, Polish – to brighten your smile, Sensitive – for gentle yet effective cleaning of sensitive gums, and White – the ideal mode to remove surface stains.

## QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

## Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

## Safe and gentle



You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns, and veneers, and helps prevent cavities and improve gum health.

## Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.



[asimpleswitch.com](https://www.asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Modes

**Clean:** For exceptional everyday clean

**White:** Removes surface stains

**Gum Care:** Gently massages gums

**Sensitive:** Gentle teeth and gum cleaning

**Polish:** Brighten and polish your teeth

## Items included

**Handle:** 1 DiamondClean

**Brush heads:** 1 W2 Optimal White, 1 C3

Premium Plaque Defense

**Glass charger:** 1

**Charger base:** 1

**Charging travel case:** 1

**Tongue brush:** 1

## Design and finishing

**Color:** Black

## Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Performance:** Removes up to 10x more plaque\*

**Whitening benefits:** Whitens teeth 2X better\*

**Timer:** Quadpacer and SmartTimer

## Ease of use

**Battery indicator:** Illuminated icon indicates battery life

**Handle:** Slim ergonomic design

**Display:** Illuminated display

**Brush head system:** Easy click-on brush heads

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** Up to 3 weeks\*

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

