

PHILIPS

sonicare

Sonic electric
toothbrush

Elite



HX9552



Custom Care system

4 ways of customized cleaning in one product

No two people are the same – especially when it comes to oral care. So why not use a power toothbrush that can adapt to your different needs and still deliver superior results.

2 brush head sizes and 2 cleaning speeds

- Standard for a complete clean, compact for precision clean
- Thorough everyday clean, massages and stimulates gums

Customizable cleaning

- Quadpacer 30-second interval timer encourages even brushing
- Easy-start to acclimate to the Sonicare experience

Superior cleaning technology

- Dynamic fluid action reaches deep between teeth
- Removes plaque in hard to reach areas
- Clinically proven safe and gentle
- Healthier gums

Highlights

Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

Healthier gums

Clinically proven to reduce gingivitis and help shrink periodontal pockets.

Specifications

Additional features

Exclusive patented Sonic technology:

Dynamic fluid action helps clean between teeth

Smartimer: Dentists' recommend a 2 minute brushing time

2 contoured brush head sizes: Standard for complete clean, Compact for precision

2 speeds: Choose the speed most comfortable for you

Easy-Start: Gently increases power over 14 uses

Programmable Quadpacer: 30-second interval timer for an even clean

Deluxe recharge gauge: Green tracks battery charge, yellow for recharge

Ergonomic tapered handle: Non-slip grip designed to improve handling

Technical specifications

Operating time (full to empty): 40

Charging time: 24 hour(s)

Battery rechargeable: NiCd

Voltage: 220-240

Power consumption: 3 W

Speed (top): Up to 31,000 movements per minute

Accessories

HX7001: Standard Elite single pack brush head

HX7002: Standard Elite double pack brush head

HX7011: Mini Elite single pack brush head

HX7012: Mini Elite double pack brush head

Weight and dimensions

F-box dimensions: 17.2 x 9.5 x 23.5 in

F-box weight: 680 lb

A-box weight: 4.08 lb

A-box dimensions: 35.1 x 29.7 x 24.1 in

A-box number of pieces: 6 pcs

Cleaning performance

Speed: Up to 31000 brush movement/min

