

# PHILIPS

## sonicare

Sonic electric  
toothbrush

Elite

HX9882



## Custom Care system

### 4 ways of customized cleaning in one product

No two people are the same - especially when it comes to oral care. So why not use a power toothbrush that can adapt to your different needs and still deliver superior results.

#### **2 brush head sizes and 2 cleaning speeds**

- Standard for a complete clean, compact for precision clean
- Thorough everyday clean, massages and stimulates gums

#### **Customizable cleaning**

- Quadpacer 30-second interval timer encourages even brushing
- Easy-start to acclimate to the Sonicare experience

#### **Superior cleaning technology**

- Dynamic fluid action reaches deep between teeth
- Removes plaque in hard to reach areas
- Clinically proven safe and gentle
- Healthier gums

# Highlights

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## Healthier gums

Clinically proven to reduce gingivitis and help shrink periodontal pockets.

# Specifications

## Additional features

**Smartimer:** Dentists' recommend a 2 minute brushing time

**2 contoured brush head sizes:** Standard for complete clean, Compact for precision

**2 speeds:** Choose the speed most comfortable for you

**Easy-Start:** Gently increases power over 14 uses

**Programmable Quadpacer:** 30-second interval timer for an even clean

**Deluxe recharge gauge:** Green tracks battery charge, yellow for recharge

**Ergonomic tapered handle:** Non-slip grip designed to improve handling

## Technical specifications

**Charging time:** Hours hour(s)

**Battery rechargeable:** NiCd

**Speed (top):** Up to 31,000 movements per minute

**Operating time (full to empty):** 28 brushings

**Voltage:** 220-240 V

## Accessories

**HX7001:** Standard Elite single pack brush head

**HX7002:** Standard Elite double pack brush head

**HX7011:** Mini Elite single pack brush head

**HX7012:** Mini Elite double pack brush head

## Weight and dimensions

**F-box dimensions:** 17.2 x 9.5 x 23.5 in

**F-box weight:** 0.68 lb

**A-box weight:** 4.08 lb

**A-box dimensions:** 35.1 x 29.7 x 24.1 in

**A-box number of pieces:** 6 pcs

## Cleaning performance

**Speed:** Up to 31000 brush movement/min

