

# PHILIPS

## sonicare


Sonic electric  
toothbrush

Essence+

1 mode

1 brush head

HX3211/65



## Helps improve gum health

Designed to be the perfect toothbrush to transition from your manual toothbrush. Until now, it feels like you've never brushed your teeth before.

### **Guides you to follow dentist recommendations**

- EasyStart program helps you adjust from a manual toothbrush
- Timers to help encourage thorough brushing

### **Designed to optimize performance**

- Click-on brush head system
- Battery life
- Sonicare dynamic cleaning action drives fluid between teeth

### **Lifts away stains for naturally whiter teeth**

- Helps whiten teeth

### **Removes more plaque than a manual toothbrush**

- Removes up to 3X more plaque along the gumline\*

### **Treats and helps prevent gingivitis**

- Better than a manual toothbrush for gingivitis

# Highlights

## Up to 3x better plaque removal




Removes up to 3x more plaque along the gumline than a manual toothbrush

## Prevents & reduces gingivitis




With superior plaque removal along the gum line, this toothbrush can help reduce gingivitis.

## Helps whiten teeth




For a bright shiny smile

## Quadpacer and SmarTimer



The Quadpacer interval timer indicates when it's time to move to another quadrant of the mouth, while the SmarTimer helps you reach the dentist recommended brushing time of two minutes.

## EasyStart



The EasyStart feature adjusts the power level of the toothbrush over the first 14 uses to help you get used to the Sonicare experience, streamlining the transition from a manual toothbrush

## Click-on brush head system



This brush head clicks on and off your brush handle for a secure fit and easy maintenance and cleaning. It fits Philips Sonicare toothbrush PowerUp Battery, Essence+/Elite+/CleanCare+ handles.

## Battery life



Lasts up to 10 days between charges.

## Sonic technology



Our sonic technology creates a unique dynamic cleaning action that drives fluid deep between teeth and along the gumline for a thorough-yet-gentle clean

# Specifications

## Modes

**Clean:** For exceptional everyday clean

## Items included

**Handles:** 1 Essence+

**Brush heads:** 1 ProResults gum health

**Charger:** 1

## Design and finishing

**Color:** Dark blue

## Cleaning performance

**Speed:** Up to 31000 brush movement/min

**Performance:** Removes up to 3x more plaque\*

**Health benefits:** Improves gum health

**Whitening benefits:** Helps maintain white teeth

**Timer:** Quadpacer and SmarTimer

## Ease of use

**Brush head system:** Easy click-on brush heads

**Brushing time:** Up to 10 days

**Battery indicator:** Light shows battery status

**Handle:** Ergonomic design, Rubber grip for easy handling

## Technical specifications

**Battery:** Rechargeable

**Battery type:** NiMH

**Operating time (full to empty):** Up to 10 days

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

