

PHILIPS

sonicare

Standard sonic
toothbrush heads

C3 Premium Plaque
Defence

2-pack

Standard size

Click-on

BrushSync mode pairing

Использование нефирменных
головочек для электрических
зубных щеток
Philips Sonicare может привести
к неправильной работе прибора
даже его поломке,
также быть небезопасным
для зубов и дёсен

HX9042/17



Our deepest cleaning brush head

Up to 10x more plaque removal* for great results

Get a clean as unique as your smile. Our C3 Premium Plaque Defence adapts to the contours of your teeth and gums thanks to soft, flexing sides that give you 4x more surface contact than a regular brush head** for a comfortable, deep clean.

For exceptional plaque removal

- Up to 10x more plaque removal than a manual toothbrush

Philips Sonicare's deep clean

- Up to 4x more surface contact** for an effortless deep clean

Innovative technology


- Philips Sonicare's advanced sonic technology
- Automatically selects optimal mode for ultimate results***

Designed to optimize performance

- Always know when to replace. Always get an effective clean.
- Click-on design for simple brush head placement
- Tested to meet your oral health needs


Highlights

Say goodbye to plaque



And thanks to its flexible design, it removes up to 10x more plaque from hard-to-reach areas* for a truly deep clean along the gum line and between teeth.

Adaptive cleaning technology



You get a personalized clean every time you brush with our adaptive cleaning technology. Soft rubber sides flex to let Premium Plaque Defence adapt to the shape of your teeth and gums, absorbing any excess brushing pressure and enhancing our sonic cleaning power. Bristles can then adjust to the shape of your teeth and gums so you get up to 4x more surface contact than a regular brush head** for deeper cleaning in hard-to-reach areas.

BrushSync™ mode-pairing



You'll always get the best possible clean with our BrushSync™ mode-pairing feature*. The Philips Sonicare C3 Premium Plaque Defence brush head syncs with your BrushSync™ enabled Philips Sonicare toothbrush handle****, selecting the optimal brushing mode and intensity level for exceptional cleaning. All you need to do is start brushing.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Brush head replacement cues



Brush heads become less effective after 3 months of use, but with BrushSync™ you'll be reminded before this happens. Your smart

toothbrush will track how often and how hard you brush, and will notify you when it's time for a replacement. Don't have a smart Philips Sonicare toothbrush? Just monitor your blue replacement bristles, and when they turn white you'll know it's time for a fresh brush head.

Brushing with peace of mind



All Philips Sonicare brush heads are safe and gentle on teeth and gums. During tests, each brush head is put through its paces to ensure durability and an exceptional performance every time you brush.

Easy to click-on



Your C3 Premium Plaque Defence brush head fits perfectly with any Philips Sonicare toothbrush handle, except PowerUp Battery and Essence. Simply click on and off for easy replacement and cleaning.

Specifications

Items included

Brush heads: 2 C3 Premium Plaque Defence

Design and finishing

Smart Brush Head Recognition

Bristle stiffness feel: Medium soft

Color: White

Material brush head: Soft, flexible rubber sides

Reminder bristles: Blue bristle color fade away

Size: Standard

Compatibility

Brush head system: Click-on

Suitable for these models: 2 Series plaque control, 2 Series plaque defense, 3 Series gum health, DiamondClean, DiamondClean Smart, EasyClean, Essence+, FlexCare, FlexCare Platinum, FlexCare Platinum Connected, FlexCare+, for Kids, HealthyWhite, HealthyWhite+, PowerUp, ProtectiveClean

Quality and performance

Tested: for optimal usage

Replacement: Every 3 months

Health benefits

Plaque removal: Removes 10x more plaque*

