

Philips
Wake-up Light

Brown/black

HF3462



Start your day in a pleasant way

The light that wakes you naturally

The Wake-up Light wakes you easily and naturally, leaving you with as much energy as on a summer's morning. Drawing on Philips' 100 years of lighting expertise, the Wake-up Light offers medically proven benefits.

Light wakes you naturally

- Gradual light positively affects energy hormones

Light intensity

- Up to 400 Lux for both natural awakening and easy reading

Natural wake-up sounds

- Three natural sounds or friendly beep

Wake up to your favourite radio station

- Digital FM radio

Philips — the light expert

- 100 years of lighting expertise

PHILIPS
sense and simplicity

Highlights

Light wakes you naturally



Light gradually increases over 30 minutes before your set wake up time. Light falls on your eyes and positively affects your energy hormones, preparing your body to wake up. This makes waking up a more pleasant feeling.

400 Lux light intensity



The Wake-up Light (HF 3461) offers up to 400 Lux light intensity. It has been clinically proven that you need up to 250 Lux to affect energy hormones in such a way that you wake up full of energy. For easy reading, you need a light level of 300 to 400 Lux.

Three natural sounds or friendly beep



At your set wake-up time, you will start hearing a natural sound or, if you prefer, a friendly beep. It starts out quietly and takes a minute and a half to reach the volume level you selected. There are three natural sounds you can choose from: morning birds in the forest, waves at the beach, or a frog and a duck in a small pond. The gradually increasing sound will

make the wake-up process nice and gentle, ensuring you're never rudely awakened. When the sound reaches your selected volume level, getting out of bed is easy.

Digital FM radio



The built-in digital FM radio lets you wake up to your favourite songs or news station every day. Just select your favourite radio station and at your set wake-up time it will start to play, at first gently, gradually increasing over ninety seconds to the volume level you selected.

Philips — the light expert

Philips is a leading company worldwide in the field of lighting technology. With more than 100 years' experience in this field, Philips is a true expert on light and its effects on the human body. The Wake-up Light combines all this expertise in a meaningful, technologically advanced and easy-to-use product that will improve the wellbeing of many.

Supported by the Medical Device Directive (MDD)



This device has proven health benefits, as certified by the Medical Device Directive (MDD).

Specifications

Technical specifications

- Power: 100 W
- Voltage: 230 V
- Frequency: 50 Hz
- Type of lamps: 1 Philips incandescent lamp
- Lifetime of lamps: 1000 hours (easy to replace)
- Insulation: class II (double isolation)
- Cord length: 3 m

Weight and dimensions

- Product dimensions: 13.5 x 20.0 x 30.0 (DxWxH) cm
- Product weight: 1.5 kg
- F-box dimensions: 16.0 x 21.0 x 32.0 (DxWxH) cm
- F-box weight: 2.1 kg
- No. of F-boxes in A-box: 3 pcs
- A-box dimensions: 22.5 x 49.5 x 34.1 (DxWxH) cm
- A-box weight: 7.8 kg
- Qnt. on Euro pallet: 126 pcs

Logistic data

- CTV code: 884346220000
- Country of origin: China

Safety

- UV-free: no UV radiation
- Cenelec certified

Easy to use

- Quick on/off button
- Snooze button for sound and light: 9 minute(s)
- Quick check function (90 sec.): Check next morning's light/sound level in 90 sec.
- Compact size: Limited bedside table space required (13.5x20 cm)
- Electricity back-up function: 5 minute(s)
- Anti-slip rubber feet

Medical appliance

- Proven health benefits: supported by Medical Device Directive certificate

Light

- Wake-up Light bulb with natural daylight spectrum: 100 W
- Sunrise simulating process (0-selected intensity): 30 minute(s)
- Light for easy reading: 400 Lux at 40–50 cm
- Light intensity dimmer: 0–400 Lux at 40–50 cm
- Light on/off function

Sound

- Wake-up sound: Three natural sounds (morning birds, seashore, pond), Friendly wake-up beep, Your favourite digital FM radio station
- Sound off function: Option to mute sound immediately
- Premium sound box: 2 W
- Digital FM radio

Display

- Easy-to-read LCD display: 50x50 mm
- Display brightness control: Three settings (low, medium, high)



Issue date 2009-05-26

Version: 3.0.1

12 NC: 0000 000 00000

© 2009 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

www.philips.com