

PHILIPS

sonicare

2 Series plaque control

1 mode

2 brush heads

2 handles



HX6212/19



Removes up to 2X more plaque* *than a manual toothbrush

Offers you superior plaque removal at the push of a button

Provides a superior clean

- Removes up to 2x more plaque than a manual toothbrush
- Better plaque removal to help reduce cavities

Designed around you

- Easy start program builds your Philips Sonicare routine
- Timers encourage thorough 2 minutes of brushing

Optimize your brushing

- Angled brush head for better access to back teeth
- Philips Sonicare's advanced sonic technology

Proven to improve oral health

- Safe and gentle for your oral health needs



asimpleswitch.com

Highlights

Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

Helps reduce cavities



Provides better plaque removal and access to back teeth; brushing twice a day every day with this brush can help reduce cavities

The only true sonic technology




Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Better access to back teeth




This plaque-removing toothbrush has a slim, angled brush head which provides better access to back teeth and hard-to-reach places.

QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

Safe and gentle



Whatever your specific needs, you can be sure of a safe brushing experience with the 2 Series plaque control. Our sonic technology is suitable for use with braces, fillings, crowns and veneers, and can even be used in the daily treatment of gum disease.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Modes

Clean: For exceptional everyday clean

Items included

Handle: 2 2-series plaque control

Brush heads: 2 C1 SimplyClean standard

Charger: 1

Design and finishing

Color: Cambridge Blue Green, Black

Cleaning performance

Speed: Up to 31000 brush movement/min

Performance: Up to 2X more plaque removal*

Health benefits: Helps reduce cavities

Timer: Quadpacer and SmarTimer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 2 weeks

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Technical specifications

Battery: Rechargeable

Battery type: NiMH

Operating time (full to empty): Up to 2 weeks

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

