

PHILIPS

sonicare

Sonic electric
toothbrush

ProtectiveClean
4300

Built-in pressure sensor

1 cleaning mode

1 x BrushSync feature



HX6888/98



Healthier gums. Gently does it.

Improve gum health up to 100% more than a manual

Feel the difference of a gentle clean with our Pressure Sensor while improving your gum health up to 100% more vs a manual toothbrush

Improves patients' oral health

- Improves gum health up to 100%*
- Fillings, crowns, or sensitive areas? Not to worry.

Innovative technology

- A gentle nudge to go nice and easy
- Patients gets the most from their brush heads
- Connects smart brush handle and smart brush heads

Proven to improve oral health

- Philips Sonicare's advanced sonic technology

Designed around your patient


- A great way to get used to your new toothbrush
- Encouragement to brush thoroughly



asimpleswitch.com

Highlights

Improves gum health



Densely-packed, high-quality bristles give you an extra gentle brushing to remove plaque along the gumline to improve gum health up to 100% more than a manual toothbrush. And, the specially curved power tip makes reaching the teeth at the back of your mouth a breeze.

Too much pressure?



Patients may not notice if they are brushing too hard, but their toothbrush will. If patients start taking it out on their gums, the toothbrush will make a pulsing sound as a reminder to ease off on the pressure.

Safe and gentle




You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns, and veneers, and helps prevent cavities and improve gum health.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Smart replacement reminders




All brush heads wear out over time. But we've got that covered too: our smart brush-recognition technology tracks how long patients have been using their brush head for, and how hard they've been brushing. When it's time to replace it, a light on the handle lets them know. That way, they can be sure that their brush head is consistently doing a stellar job.

BrushSync technology




A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

Helpful timers



Our QuadPacer lets patients know when they've spent just the right amount of time cleaning each part of their mouth, while our Smartimer tells them when they've brushed for the recommended two minutes.

Easy-start program



Starting something new can take getting used to. That's why our easy-start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Modes

Clean: For exceptional everyday clean

Items included

Handle: 1 ProtectiveClean

Brush heads: 1 G3 Premium Gum Care, 1 InterCare standard, 1 W2c Optimal White compact

Charger: 1

Design and finishing

Color: Mid-Blue

Cleaning performance

Speed: Up to 62000 brush movement/min

Health benefits: Improves gum health in only two weeks, Helps reduce cavities

Pressure feedback: Vibrates handle to alert user

Timer: Quadpacer and SmarTimer

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

Ease of use

Handle compatibility: Easy click-on brush heads

Replacement reminder: To always ensure best results, reminder icon lights up

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Brushing time: Up to 2 weeks

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 2 weeks

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

