

PHILIPS

sonicare

Sonic electric
toothbrush

DiamondClean

5 modes

2 brush heads

Glass charger



HX9342/04



Whiter, healthier teeth for life

Best Sonicare whitening electric toothbrush

Sonicare's best whitening in our most elegant Philips Sonicare electric toothbrush. Switch to Sonicare.

Proven to improve oral health

- Removes up to 7x more plaque than a manual toothbrush
- Improves gum health in only two weeks
- Whitens teeth more than 2x better than a manual toothbrush
- Sensitive Mode: Gentle teeth and gum cleaning
- Philips Sonicare toothbrush helps whiten teeth

Provides a superior clean

- DiamondClean toothbrush head for Sonicare's best whitening
- Philips Sonicare's advanced sonic technology

Designed around you


- Dual Charging System: Charging glass and travel case.

Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

Highlights

Removes up to 7x more plaque



Removes up to 7 times more plaque from hard to reach places than a manual toothbrush.

Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

Whitens teeth 2X better



This Philips Sonicare toothbrush whitens teeth more than 2x better than a manual toothbrush in just 1 week. Removes up to 100% more stains for whiter teeth in just 1 week.

DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond shaped bristles to effectively yet gently remove surface stains and scrape away plaque. Removes stain build up from coffee, tea, red wine, cigarettes and tobacco in just 1 week.

The only true sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Specifications

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Removes up to 100% more plaque

Health benefits: Improves gum health in only two weeks

Whitening benefits: Whitens teeth 2X better*

Timer: Quadpacer and Smartimer

Ease of use

Brush head system: Easy snap-on brush heads for optimal hygiene

Brushing time: Up to 3 weeks**

Battery indicator: Illuminated icon indicates battery life

Handle: Slim ergonomic design

Display: Illuminated display

Items included

Brush heads: 2 DiamondClean

Technical specifications

Battery: Rechargeable

Power consumption: 2.0 W

Battery type: Lithium ION

Design and finishing

Color: Ceramic white

Service

Warranty: 2-year limited warranty

Power

Voltage: 110-220 V

Modes

Clean: For exceptional everyday clean

Gum Care: Gently massages gums

Polish: Brighten and polish your teeth

Sensitive: Gentle teeth and gum cleaning

White: Removes surface stains

