

PHILIPS

sonicare

Interdental -
Rechargeable

AirFloss

with Fill & Charge Station

w/ProtectiveClean toothbrush

1 nozzle

1 brush head



HX8272/20



Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard to reach areas.

Improves oral health

- Up to 99.9% plaque removal
- Helps prevent cavities between teeth
- Cleans your entire mouth in just 30 seconds

The easiest way to effectively clean between teeth

- Fresh and clean with mouthwash
- An easy way to start a healthy habit

Superior technology


- Air and micro-droplet technology

Optimize your brushing

- Choose from three modes
- Improves gum health up to 100%*
- Connects smart brush handle and smart brush heads
- Timers encourage thorough 2 minutes of brushing


Highlights

Up to 99.9% Plaque Removal




AirFloss uses microdroplets of air and water or mouthwash to remove up to 99% more plaque in-between teeth*

Helps prevent cavities




By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

Fresh and clean with mouthwash




Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

30 Second Clean



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from space to space. It will automatically burst every second.

Three modes




This toothbrush allows you to customize your brushing to your own needs, with a choice of three modes. Clean mode is standard for superior cleaning. White is the ideal mode for removing surface stains. And Gum Care mode adds an extra minute of reduced-power brushing, so you can gently massage your gums.

BrushSync technology



A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

Improves gum health



Densely-packed, high-quality bristles give you an extra gentle brushing to remove plaque along the gumline to improve gum health up to 100% more than a manual toothbrush. And, the specially curved power tip makes reaching the teeth at the back of your mouth a breeze.

Specifications

Ease of use

Handle: Slim ergonomic design

Battery indicator: Flashing indicates low battery

Battery Life: Up to 2 weeks, Up to 3 wks when fully charged

Nozzle attachment: Easily snaps on and off

Technical specifications

Battery: Rechargeable

Operating time: Replace nozzle every 6 months

Battery type: Lithium ION

Charge time: 24 hours to full charge

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty

Items included

AirFloss Handle: 1

Fill & Charge station: 1

AirFloss Nozzle: 1

Handle: 1 ProtectiveClean

Brush heads: 1 G2 Optimal Gum Care

Cleaning performance

Health benefits: Helps improve gum health

Auto-burst feature: Cleans entire mouth in 30 sec

Speed: Up to 62000 brush movement/min

Design and finishing

Color: White and transparent green, White and Mint

Modes

Clean: For exceptional everyday clean

Gum Care: Gently massages gums

White: Removes surface stains

