

# PHILIPS

## sonicare

Sonic electric  
toothbrush

ProtectiveClean  
4300

Built-in pressure sensor

1 cleaning mode

1 x BrushSync feature



HX6809/16



## Healthier gums. Gently does it.

Improve gum health up to 100% more than a manual

Feel the difference of a gentle clean with our Pressure Sensor while improving your gum health up to 100% more vs a manual toothbrush

### Proven to improve oral health

- Improves gum health by up to 100%\*
- Safe & gentle on sensitive areas, orthodontics & dental work

### Innovative technology

- Lets you know if you're pressing too hard
- Always know when to replace your brush heads
- Philips Sonicare's advanced sonic technology
- Connects smart brush handle and smart brush heads

### Designed around you

- Easy start program builds your Philips Sonicare routine
- Encouragement to brush thoroughly


# Highlights

## Improves gum health



Densely-packed, high-quality bristles give you an extra gentle brushing to remove plaque along the gumline to improve gum health up to 100% more than a manual toothbrush. And, the specially curved power tip makes reaching the teeth at the back of your mouth a breeze.

## Monitor your brushing pressure



Brushing too hard can damage your teeth and gums. To prevent this, your Philips Sonicare ProtectiveClean emits a gentle pulsing sound to remind you to ease off on the pressure.

## Safe and gentle



You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns, and veneers, and helps prevent cavities and improve gum health.

## The only true sonic technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

## BrushSync replacement reminder




All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. When it's time to replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

## BrushSync technology




A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

## Helpful timers



Need an electric toothbrush with a timer? Our QuadPacer lets you know when you've spent just the right amount of time cleaning each part of your mouth, while our Smartimer tells you when you've brushed for the recommended two minutes.

## Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

# Specifications

## Modes

**Clean:** For exceptional everyday clean

**2 intensities:** Low, High

## Items included

**Handle:** 1 ProtectiveClean

**Brush heads:** 1 G2 Optimal Gum Care

**Charger:** 1

## Design and finishing

**Color:** White and Light Blue

## Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Health benefits:** Improves gum health in only two weeks, Helps reduce cavities

**Pressure feedback:** Vibrates handle to alert user

**Timer:** Quadpacer and SmarTimer

## Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**BrushSync Replacement Reminder:** Always know when to, replace brush heads

## Ease of use

**Handle compatibility:** Easy click-on brush heads

**Battery indicator:** Light shows battery status

**Handle:** Slim ergonomic design

**Brushing time:** Up to 2 weeks

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** Up to 2 weeks

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

