

# PHILIPS

## sonicare

2 Series plaque control

EssentialClean

1 mode

1 brush head

HX6231/72



## Removes up to 7X more plaque\*

\*than a manual toothbrush

Offers you superior plaque removal at the push of a button

#### **Proven to improve oral health**

- Up to 7x more plaque removal than a manual toothbrush
- Safe and gentle for your oral health needs

#### **Innovative technology**

- Philips Sonicare's advanced sonic technology

#### **Designed around you**

- Easy start program builds your Philips Sonicare routine
- Timers encourage thorough 2 minutes of brushing

#### **Provides a superior clean**

- Better plaque removal to help reduce cavities
- Angled neck to help tackle cleaning in tough areas

# Highlights

## Say goodbye to plaque



Densely-packed, high-quality bristles remove up to 7x more plaque than a manual toothbrush.

## Helps reduce cavities



Provides better plaque removal and access to back teeth; brushing twice a day every day with this brush can help reduce cavities

## The only true sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Safe and gentle



Whatever your specific needs, you can be sure of a safe brushing experience with the 2 Series plaque control. Our sonic technology is suitable for use with braces, fillings, crowns and veneers, and can even be used in the daily treatment of gum disease.

## Reach those tricky spots




The unique shape of your 2 Series toothbrush handle combined with our brush heads means that even tricky spots like back teeth get a thorough clean.

## Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

## QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

# Specifications

## Modes

**Clean:** For exceptional everyday clean

## Items included

**Handle:** 1 2-Series

**Brush heads:** 1 W DiamondClean

**Travel case:** 1

**Charger:** 1

## Design and finishing

**Color:** Ultra Coral

## Cleaning performance

**Performance:** Removes up to 7x more plaque\*

**Health benefits:** Helps reduce cavities

**Timer:** Quadpacer and SmartTimer

## Ease of use

**Battery indicator:** Light shows battery status

**Handle:** Slim ergonomic design

**Brush head system:** Easy click-on brush heads

**Brushing time:** Up to 2 weeks

## Technical specifications

**Battery:** Rechargeable

**Battery type:** NiMH

**Operating time (full to empty):** Up to 2 weeks

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

