

PHILIPS

sonicare

Sonic electric
toothbrush

HealthyWhite

2 modes

3 brush heads



HX6713/11



Whiter, healthier teeth

The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Philips Sonicare HealthyWhite electric toothbrush is proven to remove everyday stains in just two weeks, with the use of the Clean & White Mode.

Proven to improve oral health

- Philips Sonicare toothbrush helps whiten teeth

Provides a superior clean

- Removes up to 7x more plaque than a manual toothbrush
- DiamondClean brush head for your best whitening
- Philips Sonicare's advanced sonic technology

Offers a customized brushing experience

- Clean and White Mode: proven to remove stains
- Easy-start program to acclimate to the Sonicare experience

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing



asimpleswitch.com

Highlights

Helps whiten teeth




This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*

DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond shaped bristles that are medium stiff to effectively yet gently scrape away plaque. This Philips Sonicare electric toothbrush provides an exceptional cleaning and whiter teeth compared to a manual toothbrush.

Removes up to 7x more plaque




Removes up to 100% more plaque from hard to reach places than a manual toothbrush.

The only true sonic technology




Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

SmarTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Modes

Clean: For exceptional everyday clean

Clean & White: 2 min Clean + 30 sec White

Items included

Handles: 1 HealthyWhite

Brush heads: 1 DiamondClean standard, 1 ProResults gum health, 1 ProResults plaque control

Charger: 1

Travel case: 1

Design and finishing

Color: Glacier Green

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Removes up to 7x more plaque*

Health benefits: Helps improve gum health

Whitening benefits: Helps naturally remove stains, Whitens teeth up to 2 shades

Timer: Quadpacer and SmarTimer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 3 weeks**

Battery type: Lithium ION

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

