

# PHILIPS

## sonicare

Sonic electric  
toothbrush

ProtectiveClean  
5100

Built-in pressure sensor

3 modes

1 x BrushSync feature

Travel case



HX6859/68

## Healthier gums. Gently does it.

Improve gum health up to 100% more than a manual

Feel the difference of a gentle clean with our Pressure Sensor while improving your gum health up to 100% more vs a manual toothbrush

### Proven to improve oral health

- Safe & gentle on sensitive areas, orthodontics & dental work
- Whitens teeth in just one week

### Optimize your brushing

- Choose from three modes

### Innovative technology

- Lets you know if you're pressing too hard
- Always know when to replace your brush heads
- Philips Sonicare's advanced sonic technology
- Connects smart brush handle and smart brush heads

### Designed around you

- Encouragement to brush thoroughly
- Makes traveling easier for you


# Highlights

## Naturally whiter teeth



Click on the W2 Optimal White brush head to remove surface stains and reveal a whiter smile. With its densely packed central stainremoval bristles, it's clinically proven to whiten teeth in just one week.

## Three modes



This toothbrush allows you to customize your brushing to your own needs, with a choice of three modes. Clean mode is standard for superior cleaning. White is the ideal mode for removing surface stains. And Gum Care mode adds an extra minute of reduced-power brushing, so you can gently massage your gums.

## Safe and gentle




You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns, and veneers, and helps prevent cavities and improve gum health.

## Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

## Monitor your brushing pressure



Brushing too hard can damage your teeth and gums. To prevent this, your Philips Sonicare ProtectiveClean emits a gentle pulsing sound to remind you to ease off on the pressure.

## BrushSync technology



A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

## BrushSync replacement reminder



All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. When it's time to replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

## Helpful timers



Need an electric toothbrush with a timer? Our QuadPacer lets you know when you've spent just the right amount of time cleaning each part of your mouth, while our Smartimer tells you when you've brushed for the recommended two minutes.

## Handy travel case and charger



Our premium travel case lets you store your toothbrush hygienically, while our compact charging base keeps you topped up when you're on the go. You do get two weeks of regular use from a single full charge, but the charger's just the ticket for longer trips.

# Specifications

## Modes

**Clean:** For exceptional everyday clean

**Gum Care:** Gently massages gums

**White:** Removes surface stains

## Items included

**Handle:** 1 ProtectiveClean

**Brush heads:** 1 W2 Optimal White

**Travel case:** 1

**Charger:** 1

**UV sanitizer:** 1

## Design and finishing

**Color:** White and Light Blue

## Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Health benefits:** Improves gum health in only two weeks, Helps reduce cavities

**Pressure feedback:** Vibrates handle to alert user

**Timer:** Quadpacer and SmartTimer

## Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**BrushSync Replacement Reminder:** Always know when to, replace brush heads

## Ease of use

**Handle compatibility:** Easy click-on brush heads

**Replacement reminder:** To always ensure best

results, reminder icon lights up

**Battery indicator:** Light shows battery status

**Handle:** Slim ergonomic design

**Brushing time:** Up to 2 weeks

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** Up to 2 weeks

## Power

**Voltage:** 100-240 V, 50-60 Hz

## Service

**Warranty:** 2-year limited warranty

