

PHILIPS

sonicare

Sonic electric
toothbrush

For Kids

2 modes

1 brush head

8 Stickers



HX6311/07



Sonic power maximises your child's routine Electric toothbrush for kids

Philips Sonicare For Kids is an electric power toothbrush for kids aged 4 and older. It offers maximum plaque removal, sonic technology, customisable stickers and educational tools to help make proper brushing fun for a lifetime.

Designed to help build proper brushing habits

- KidTimer helps to increase brushing time
- 2 Kid-friendly power modes deliver a gentle, effective clean
- 2 brush head sizes available
- KidPacer delivers fun and engaging musical tones
- Fun interchangeable stickers for customisation
- Multi grip design for parents and kids
- Preferred by dental professionals for their own children
- Anti-roll shape

Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Removes more plaque than a children's manual toothbrush

Safe for growing smiles

- Rubberised brush head is designed to protect young teeth


Highlights

Age-appropriate brush heads



This Philips Sonicare electric toothbrush has two brush head sizes that are specifically designed to gently clean and protect teeth as they grow

KidPacer



Fun sounds encourage kids to brush their entire mouth effectively for the dentist-recommended 2 minutes

KidTimer



To help kids ease into a proper routine, this electric toothbrush slowly increases the brushing time over 90 days until it reaches the dentist-recommended 2 minutes to help instil healthy habits naturally.

Power modes




With two kid-friendly power modes, this power toothbrush provides proper cleaning for different ages; a low mode for younger children and a high mode for older kids.

Stickers



Allows kids to make each brush uniquely their own

Maximum plaque removal



Philips Sonicare For Kids removes up to 75% more plaque than manual toothbrushes*

Multi grip design




Ergonomically designed to help kids brush properly on their own

Sonic technology



Features a unique dynamic action that gently and effectively reaches deep between teeth and along the gumline.

Preferred by the dental profession



91% of dental professional parents prefer Philips Sonicare For Kids for their own children**

Rubberised brush head



The rubberised brush head of this electric toothbrush is designed to help protect growing teeth

Specifications

Modes

Power modes: 2

Items included

Handles: 1 Sonicare for Kids

Brush heads: 1 Sonicare for Kids standard

Charger: 1

Stickers: 8 customisation stickers, 2 Bonus stickers

Design and finishing

Colour: Aqua

Cleaning performance

Speed: Up to 31,000 brush strokes/min

Performance: 75% more effective*

Health benefits: For healthy oral care habits

Timer: KidTimer and Quadpacer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Light shows battery status

Handle: Slim ergonomic design, Rubber grip for easy handling

Display: Illuminated display

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 3 weeks**

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty



* than a manual toothbrush

* in hard-to-reach areas

** Philips Sonicare in a home survey of US dental professionals with children aged 4-10

*** based on two, two-minute periods of brushing per day