

# PHILIPS

## sonicare

Compact sonic  
toothbrush heads

ProResults


2-pack

Compact size

Click-on

All-around cleaning

HX6022/08



## A focused clean with Philips Sonicare ProResults

### Exceptional performance for healthy teeth and gums

With contoured bristles to follow the shape of teeth and gums, Philips Sonicare ProResults brush heads use precision engineering and the power of sonic motion to help patients brush in hard-to-reach areas for a complete, gentle clean.

#### **Designed to optimize performance**

- Superior bristle design
- Reminder bristles ensure your most effective clean

#### **Proven to improve oral health**

- Part of a better oral healthcare routine

#### **Provides a superior clean**

- Engineered to maximize sonic motion

#### **Designed around you**

- Works with any Philips Sonicare click-on toothbrush


# Highlights

## Superior bristle design



Philips Sonicare ProResults brush heads use contoured and angled bristles to follow the shape of teeth and maintain surface contact, maximizing patients' ability to remove plaque in places they might normally miss with a manual toothbrush. With color-fading reminder bristles to indicate effectiveness, patients will always know when it's time for a new brush head.

## Reminder bristles



At first glance it may not be obvious, but brush heads lose stiffness and gradually wear down over months of normal use. Our blue reminder bristles fade to white and help you recognize when it's time for a replacement. For optimal results, replace your brush head every three months.

## Click-on design



Philips Sonicare brush heads click on and off the brush handle for a secure fit and easy maintenance and cleaning. Fits all Philips Sonicare toothbrush handles except: PowerUp and Essence.

## Better oral healthcare



Like all authentic Philips Sonicare branded brush heads, this brush head is safe on teeth and gums. Each brush head has been quality tested for exceptional performance and durability.

## Sonic technology



Philips Sonicare brush heads are vital to our core technology of high-frequency, high-amplitude brush movements performing over 31,000 brush strokes per minute. This dynamic cleaning action drives fluid deep between teeth and along the gumline, for a thorough yet gentle cleaning.

# Specifications

## Items included

Brush heads: 2 ProResults compact

## Compatibility

Brush head system: Click-on

Suitable for these models: 2 Series plaque control, 3 Series gum health, DiamondClean, EasyClean, FlexCare, FlexCare Platinum, FlexCare+, for Kids, HealthyWhite, HealthyWhite+

## Design and finishing

Bristle stiffness feel: Medium

Color: White

Reminder bristles: Blue bristle color fade away

Size: Compact

## Health benefits

Plaque removal: Helps remove plaque

## Quality and performance

Replacement: Every 3 months

Tested: for optimal usage

