

# PHILIPS

## sonicare

Sonic electric  
toothbrush

HealthyWhite

2 modes

2 brush heads



HX6762/20



## Whiter, healthier teeth

Removes up to 100% more stains for whiter teeth\*

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Sonicare HealthyWhite electric toothbrush is proven to remove everyday stains in just two weeks, with the regular use of the Clean & White Mode.

### Provides a superior clean

- DiamondClean toothbrush head for Sonicare's best whitening
- Removes up to 7x more plaque than a manual toothbrush
- Philips Sonicare's advanced sonic technology
- Angled brush head neck for better reach to back teeth

### Proven to improve oral health

- Naturally whiter teeth with patented sonic technology
- Clinically proven safe and gentle

### Offers a customized brushing experience

- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2 minute mode
- Easy-start program to acclimate to the Sonicare experience

### Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

# Highlights

## Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

## Clean and White Mode




2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.\*

## Clean Mode

For exceptional everyday clean. Helps maintain whitening treatments.

## DiamondClean brush head




DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond shaped bristles to effectively yet gently remove surface stains and scrape away plaque. Removes stain build up from coffee, tea, red wine, cigarettes and tobacco in just 1 week.

## The only true sonic technology




Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Removes up to 7x more plaque



Removes up to 7 times more plaque from hard to reach places than a manual toothbrush.

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Angled brush head neck




The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## SmartTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

# Specifications

## Modes

**Clean:** For exceptional everyday clean

**Clean & White:** 2 min Clean + 30 sec White

## Items included

**Handles:** 1 HealthyWhite

**Brush heads:** 2 DiamondClean standard

**Charger:** 1

## Design and finishing

**Color:** Pink

## Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Performance:** Removes up to 7x more plaque\*,  
Up to 100% more stain removal\*

**Health benefits:** Helps improve gum health

**Whitening benefits:** Whitens teeth up to 2  
shades

**Timer:** Quadpacer and SmarTimer

## Ease of use

**Brush head system:** Easy click-on brush heads

**Brushing time:** Up to 3 weeks\*\*

**Battery indicator:** Light shows battery status

**Handle:** Slim ergonomic design

## Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** Up to 3  
weeks\*\*

**Battery type:** Lithium ION

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

