

# PHILIPS

## sonicare

Sonic electric  
toothbrush

FlexCare+

5 modes

1 brush head

HX6921/43



## Complete gum care

Toothbrush for gentle cleaning of teeth and gums

Combining our sonic technology and 5 brushing modes, FlexCare+ electric toothbrush is your solution for complete gum care. Now you can improve your gum health after just 2 weeks of use.

### Proven to improve oral health

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology
- Improves gum health in only two weeks

### Offers a customized brushing experience

- Five different modes for a better cleaning experience
- Easy-start program to acclimate to the Sonicare experience

### Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

### Provides a superior clean


- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth



[asimpleswitch.com](http://asimpleswitch.com)


# Highlights

## Five different cleaning modes



The five modes include a Clean mode for an outstanding clean in two minutes. A Sensitive mode for gentle yet thorough cleaning of the gums. A Refresh mode for those quick one minute touch-ups. A Massage mode to stimulate the gum and a Gum Care mode that helps improve gum health.

## Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck




The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## SmarTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Modes

**Clean:** For exceptional everyday clean

**Gum Care:** Gently massages gums

**Refresh:** 1 min mode for quick touch up

**Sensitive:** Gentle teeth and gum cleaning

**Massage:** Invigorating massage

## Items included

**Handles:** 1 FlexCare+

**Brush heads:** 1 ProResults gum health

**Charger:** 1

## Design and finishing

**Color:** Cooper frost

## Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Performance:** Removes up to 6X more plaque\*

**Health benefits:** Helps improve gum health

**Whitening benefits:** Helps naturally remove stains

**Timer:** Quadpacer and SmarTimer

## Ease of use

**Brush head system:** Easy click-on brush heads

**Brushing time:** Up to 3 weeks\*\*

**Battery indicator:** Light shows battery status

**Handle:** Slim ergonomic design

**Display:** Illuminated display

## Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** Up to 3 weeks\*\*

**Battery type:** Lithium ION

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

