

PHILIPS

sonicare

Sonic electric
toothbrush

DiamondClean

5 modes

3 brush heads

Glass charger

USB travel charger



HX9333/07



Whiter, healthier teeth for life

Best Philips Sonicare whitening sonic toothbrush

Sonicare's best whitening in our most elegant Philips Sonicare electric toothbrush. Switch to Sonicare.

Proven to improve oral health

- Up to 7x more plaque removal than a manual toothbrush
- A whiter smile in 1 week with our DiamondClean brush head*
- Improves gum health in just 2 weeks*
- Choose from 5 modes including Polish and Sensitive
- Safe and gentle for your oral health needs

Provides a superior clean

- Our unique technology brings you a powerful yet gentle clean
- Angled neck to help tackle cleaning in tough areas

Designed around you

- Timers encourage thorough 2 minutes of brushing
- Easy start program builds your Philips Sonicare routine
- USB charging travel case and charging glass



asimpleswitch.com

Specifications

Modes

Clean: For exceptional everyday clean

Gum Care: Gently massages gums

Polish: Brighten and polish your teeth

Sensitive: Gentle teeth and gum cleaning

White: Removes surface stains

Items included

Handle: 1 DiamondClean

Brush heads: 1 DiamondClean standard, 1 InterCare standard, 1 Sensitive standard

Glass charger: 1

Travel case: USB travel charger

Tongue scraper: 1

Design and finishing

Color: Ceramic white

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Removes up to 7x more plaque*

Health benefits: Improves gum health in only two weeks

Whitening benefits: Whitens teeth 2X better*

Timer: Quadpacer and SmarTimer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Illuminated icon indicates

battery life

Handle: Slim ergonomic design

Display: Illuminated display

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 3 weeks**

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

* **based on two periods of two-minute brushings per day, on clean mode