

PHILIPS

sonicare

AirFloss Pro/Ultra -
Interdental cleaner

DiamondClean

w/AirFloss Ultra cleaner

1 nozzle

1 brush head

Travel case



HX8492/01



Whiter, healthier teeth for life

Best Philips Sonicare whitening sonic toothbrush

Philips Sonicare's best whitening in our most elegant sonic electric toothbrush. Switch to Philips Sonicare.

Provides a superior clean

- Up to 10x more plaque removal than a manual toothbrush
- Philips Sonicare's advanced sonic technology
- Choose from 5 modes including Deep Clean and Sensitive
- Deep Clean mode tackles your trouble spots
- USB charging travel case and charging glass

Improves oral health

- Whitens teeth more than 2x better than a manual toothbrush
- Up to 99.9% plaque removal***

Superior technology


- Air and micro-droplet technology
- Triple burst customizable settings

The easiest way to effectively clean between teeth

- Takes just 60 seconds: Point, Press, Clean!


Highlights

Say goodbye to plaque



And thanks to its flexible design, it removes up to 10x more plaque from hard-to-reach areas* for a truly deep clean along the gum line and between teeth.

Whitens teeth 2X better



This Philips Sonicare toothbrush whitens teeth more than 2x better than a manual toothbrush in just 1 week. Removes up to 100% more stains for whiter teeth in just 1 week.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

5 modes to clean your way



With DiamondClean, you get a truly refreshing clean every day. Our 5 modes cover all your brushing needs; Clean mode – for exceptional daily cleaning, Gum Care – to gently massage gums, Deep Clean – for an invigorating deep clean, Sensitive – for gentle yet effective cleaning of sensitive gums, and White – the ideal mode to remove surface stains.

Extra deep cleaning



Our Deep Clean mode means you can really focus on specific areas, paying them extra attention. We've combined our unique brush head motion with an extended 3 minutes of brushing time so your teeth get a truly deep clean.

Design with storage in mind



Your deluxe USB travel case doubles as a charger so you stay topped up on the go. Place your toothbrush in the case and connect to your laptop or plug into a wall socket. Also includes a brush head holder for extra hygiene


when traveling. As for charging at home, our sleek charger glass fits stylishly into your bathroom and doubles as a glass for rinsing after brushing. Simply place your brush into the glass. Enjoy two weeks of regular use from a single full charge.

Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.***

Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

Specifications

Cleaning performance

Clinically proven results: As effective as floss, for gum health improvement, for plaque reduction, for reducing gum bleeding

Customizable settings: Auto-burst or manual burst, Single/double/triple burst***

For best results: Change nozzle every 6 months

Health benefits: Cleans entire mouth in 60 sec, Improve gum health in 2 weeks*, Improves gum health in only two weeks

Performance: Removes up to 99.9% plaque**

Whitening benefits: Whitens teeth 2X better*

Speed: Up to 62000 brush movement/min

Timer: Quadpacer and SmarTimer

Design and finishing

Color: White

Ease of use

Battery indicator: Handle LEDs show charge level, Illuminated icon indicates battery life

Battery Life: 1-3 weeks*****

Handle: Slim ergonomic design

Nozzle attachment: Easily snaps on and off

Brush head system: Easy click-on brush heads

Display: Illuminated display

Items included

AirFloss Pro/Ultra handle: 1

AirFloss Pro/Ultra nozzle: 1

Handle: 1 DiamondClean

Brush heads: 1 C3 Premium Plaque Control

Travel case: Premium travel case

Charger: 1

Charger base: 1

Glass charger: 1

Power

Voltage: Multi-voltage charger, 110-220 V

Service

Warranty: 2-year limited warranty

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Charge time: 24 hours to full charge

Modes

Clean: For exceptional everyday clean

Gum Care: Gently massages gums

Sensitive: Gentle teeth and gum cleaning

White: Removes surface stains

Deep Clean: For an invigorating deep clean

