

# PHILIPS

## sonicare

Sonic electric  
toothbrush

FlexCare+

5 modes

1 brush head



HX6992



## Complete Gum Care

Toothbrush for gentle cleaning of teeth and gum

Combining patented sonic technology and 5 brushing modes, FlexCare+ electric toothbrush is your solution for complete gum care. Now you can improve your gum health after just 2 weeks of use.

### Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

### Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- UV brush head Sanitizer - Cleaner brush, better oral health

### Proven to improve oral health

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology
- Improves gum health in only two weeks

### Offers a customized brushing experience

- Easy-start program to acclimate to the Sonicare experience
- Five different modes for a better cleaning experience

# Highlights

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck




The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## UV brush head Sanitizer

All-in-one sanitizing, charging and storage. The UV brush head Sanitizer helps keep your brush head free of bacteria. It kills up to 99% of bacteria and viruses.\*

## Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## Improves gum health




This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

## SmarTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Five different cleaning modes



The five modes include a Clean mode for an outstanding clean in two minutes. A Sensitive mode for gentle yet thorough cleaning of the gums. A Refresh mode for those quick one minute touch-ups. A Massage mode to stimulate the gum and a Gum Care mode that helps improve gum health.

# Specifications

## Technical specifications

**Battery:** Rechargeable

## Design and finishing

**Color:** White with sky blue grip and illuminated display

## Items included

**Charger:** 1 multi-voltage charger with UV brush head sanitizer and 1 multi-voltage travel charger

**Charger cover:** 1 charger cover with brush head holder

**Brush head color rings:** 3 color rings for easy sharing

## Ease of use

**Brush head system:** Easy snap-on brush heads for optimal hygiene

**Brushing time:** 42 2-minute brushings OR 3 weeks

**Battery indicator:** Illuminated icon indicates battery life

**Handle:** Slim ergonomic design and rubber grip for easy handling.

**Display:** Illuminated display

## Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Brush heads:** 1 ProResults standard brush head for an All-Around Clean and 1 ProResults compact brush head for a Precision Clean

**Performance:** Removes up to 83% of plaque in hard-to-reach areas

**Health benefits:** Improves gum health in only two weeks

**Whitening benefits:** Helps remove tea, coffee and other stains

**Timer:** Quadpacer and SmartTimer

**Modes:** 5 Modes: Clean, Sensitive, Massage, Refresh and Gum Care

## Service

**Warranty:** 2-year limited warranty

## Power

**Voltage:** Multi-voltage charger

