

PHILIPS

sonicare

Sonic electric
toothbrush

3 Series gum health

1 mode

3 intensity settings

2 brush heads



HX6632/24



Improves gum health up to 100%*

*more than manual brushing

Provides superior plaque removal and gum care in an easy-to-use power toothbrush to help prevent and reduce the early signs of gum disease while offering you the most comfortable brushing experience

Proven to improve oral health

- Philips Sonicare toothbrush helps whiten teeth
- Removes up to 7 x more plaque than a manual toothbrush
- Philips Sonicare's advanced sonic technology
- DiamondClean toothbrush head for Sonicare's best whitening

Offers a customised brushing experience

- Three intensity settings for a gentle experience

Recommended by dental professionals worldwide

- Proven to improve gum health in 3 ways
- Safe on orthodontics, dental work and veneers

Guides you to follow dentist recommendations

- Smartimer helps you meet recommended brushing time
- Quadpacer interval timer encourages thorough brushing



asimpleswitch.com


Highlights

Improves gum health in 3 ways



Clinically proven to improve gum health in three ways: reduces redness, swelling and bleeding (sites)

Helps whiten teeth




This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Removes up to 7 x more plaque



Removes up to 7 times more plaque from hard-to-reach places than a manual toothbrush.

Three intensity settings



With three customisable intensity settings (low, medium and high) you can select the right power level to ensure that brushing along the


gumline is comfortable. After 1.5 seconds of brushing, you will have to push the button twice to switch modes — once to pause and again to switch.

The only true sonic technology



The Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

Quadpacer



Features a Quadpacer interval timer to help encourage thorough brushing of the entire mouth. The toothbrush beeps every 30 seconds, signifying that it's time to move to another area of the mouth.

Safe and gentle




The gentle motion of this Philips Sonicare power toothbrush is safe for dental work, including orthodontics, implants and veneers, so you can brush without worrying.

Smartimer



This Philips Sonicare power toothbrush features a special Smartimer to help you meet the dentist-recommended brushing time of two minutes

DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond shaped bristles to effectively yet gently remove surface stains and scrape away plaque. Removes stain build up from coffee, tea, red wine, cigarettes and tobacco in just 1 week.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Modes

Clean: For exceptional everyday cleaning
3 intensities: Low, Medium, High

Items included

Handles: 1 3 Series gum health
Brush heads: 2 DiamondClean standard
Charger: 1
Travel case: 1

Design and finishing

Colour: Smoke White

Cleaning performance

Speed: Up to 31,000 brush movements/min
Performance: Removes up to 7x more plaque*
Health benefits: Improves gum health
Whitening benefits: Whitens teeth in 1 week
Timer: Quadpacer and Smartimer

Ease of use

Brush head system: Easy click-on brush heads
Brushing time: Up to 3 weeks**
Battery indicator: Light shows battery status
Handle: Slim ergonomic design

Technical specifications

Battery: Rechargeable
Battery type: Lithium ION
Operating time (full to empty): Up to 3 weeks**

Power

Voltage: 110–220 V

Service

Warranty: 2-year limited warranty

