

# PHILIPS

## sonicare

Sonic electric  
toothbrush

DiamondClean

5 modes

6 brush heads

Glass charger

USB travel charger



HX9336/25



## Whiter, healthier teeth for life

### Best Sonicare whitening electric toothbrush

Sonicare's best whitening in our most elegant Philips Sonicare electric toothbrush. Switch to Sonicare.

#### Proven to improve oral health

- A whiter smile in 1 week with our DiamondClean brush head\*
- Up to 7x more plaque removal than a manual toothbrush
- Improves gum health in just 2 weeks\*
- Choose from 5 modes including Polish and Sensitive
- Safe and gentle for your oral health needs

#### Provides a superior clean

- Our unique technology brings you a powerful yet gentle clean
- Angled neck to help tackle cleaning in tough areas

#### Designed around you

- Easy start program builds your Philips Sonicare routine
- USB charging travel case and charging glass
- Timers encourage thorough 2 minutes of brushing



[asimpleswitch.com](http://asimpleswitch.com)


# Highlights

## Say goodbye to plaque




Densely-packed, high-quality bristles remove up to 7x more plaque than a manual toothbrush.

## Up to 2x whiter teeth




Click on our DiamondClean brush head to gently yet effectively remove surface stains. Densely packed central stain removal bristles work hard to make your smile 2x whiter in just 7 days.\*

## Look after those gums



With optimal cleaning from your DiamondClean, gums will become healthier in 2 weeks\*. Removing up to 7x more plaque along the gum line than a manual toothbrush\*, you'll get your healthiest smile.

## 5 modes to clean your way



With DiamondClean, you get a refreshing clean every day. Our 5 modes cover all your brushing needs; Clean mode – for exceptional daily cleaning, Gum Care – to gently massage gums, Polish – to brighten your smile, Sensitive – for gentle yet effective cleaning of sensitive gums, and White – the ideal mode to remove surface stains.

## Design with storage in mind



Your deluxe USB travel case doubles as a charger so you stay topped up on the go. Place your toothbrush in the case and connect to your laptop or plug into a wall socket. Also includes a brush head holder for extra hygiene when traveling. As for charging at home, our sleek charger glass fits stylishly into your bathroom and doubles as a glass for rinsing after brushing. Simply place your brush into the glass. Enjoy two weeks of regular use from a single full charge.

## Philips Sonicare technology



Powerful sonic vibrations whip up your toothpaste into plaque fighting bubbles, and drive them deep between your teeth and along your gum line. At the same time your teeth will


experience 62,000 gentle but effective brush movements. You'll get a month's worth of manual brushing in just 2 minutes.

## QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

## Reach those tricky spots



The unique shape of your DiamondClean toothbrush handle combined with our brush heads means that even tricky spots like back teeth get a thorough clean.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Modes

**Clean:** For exceptional everyday clean

**Gum Care:** Gently massages gums

**Polish:** Brighten and polish your teeth

**Sensitive:** Gentle teeth and gum cleaning

**White:** Removes surface stains

## Items included

**Brush heads:** 1 DiamondClean standard, 1 DiamondClean compact, 1 InterCare standard, 1 Sensitive standard, 1 ProResults gum health, 1 ProResults plaque control

**Glass charger:** 1

**Travel case:** USB travel charger

**Handle:** 1 DiamondClean

## Design and finishing

**Color:** Ceramic white

## Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Performance:** Removes up to 7x more plaque\*

**Health benefits:** Improves gum health in only two weeks

**Whitening benefits:** Whitens teeth 2X better\*

**Timer:** Quadpacer and SmarTimer

## Ease of use

**Brush head system:** Easy click-on brush heads

**Brushing time:** Up to 3 weeks\*\*

**Battery indicator:** Illuminated icon indicates

battery life

**Handle:** Slim ergonomic design

**Display:** Illuminated display

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** Up to 3 weeks\*\*

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

