

# PHILIPS

## sonicare

Sonic electric  
toothbrush

3 Series gum health

1 mode

3 intensity settings

1 brush head



HX6631/13



## Improves gum health up to 100%\*

\*more than manual brushing

Provides superior plaque removal and gum care in an easy-to-use power toothbrush to help prevent and reduce the early signs of gum disease while offering you the most comfortable brushing experience

### Improves gum health more vs. a manual toothbrush

- Sonicare dynamic sonic action drives fluid between teeth
- Removes up to 6 X more plaque along the gumline
- Better than a manual toothbrush for gingivitis

### Guides you to follow dentist recommendations

- Smartimer helps you meet recommended brushing time
- Quadpacer interval timer encourages thorough brushing

### Offers a customised brushing experience

- Three intensity settings for a gentle experience

### Recommended by dental professionals worldwide

- Proven to improve gum health in 3 ways
- Safe on orthodontics, dental work and veneers



[asimpleswitch.com](http://asimpleswitch.com)


# Highlights

## Improves gum health in 3 ways




Clinically proven to improve gum health in three ways: reduces redness, swelling and bleeding (sites)

## Prevents and reduces gingivitis



With superior plaque removal along the gumline, this toothbrush can help reduce gingivitis up to 100% better than a manual toothbrush in just two weeks

## Three intensity settings



With three customisable intensity settings (low, medium and high) you can select the right power level to ensure that brushing along the gumline is comfortable. After 1.5 seconds of brushing, you will have to push the button twice to switch modes — once to pause and again to switch.

## Up to 6x better plaque removal



Removes up to 6 X more plaque along the gumline than a manual toothbrush.

## Sonic technology



Our sonic technology creates a unique dynamic cleaning action that drives fluid deep between teeth and along the gumline for a thorough-yet-gentle clean

## Safe and gentle



The gentle motion of this Philips Sonicare power toothbrush is safe for dental work, including orthodontics, implants and veneers, so you can brush without worrying.

## Smartimer



This Philips Sonicare power toothbrush features a special Smartimer to help you meet the dentist-recommended brushing time of two minutes

## Quadpacer



Features a Quadpacer interval timer to help encourage thorough brushing of the entire mouth. The toothbrush beeps every 30 seconds, signifying that it's time to move to another area of the mouth.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Modes

**Clean:** For exceptional everyday cleaning

**3 intensities:** Low, Medium, High

## Items included

**Handles:** 1 3 Series gum health

**Brush heads:** 1 ProResults gum health

**Charger:** 1

## Design and finishing

**Colour:** Arctic Mist

## Cleaning performance

**Speed:** Up to 31,000 brush movements/min

**Performance:** Removes up to 6 X more plaque\*

**Health benefits:** Improves gum health

**Timer:** Quadpacer and Smartimer

## Ease of use

**Brush head system:** Easy click-on brush heads

**Brushing time:** Up to 3 weeks\*\*

**Battery indicator:** Light shows battery status

**Handle:** Slim ergonomic design

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** Up to 3 weeks\*\*

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

