

PHILIPS

sonicare

Sonic electric
toothbrush

HealthyWhite+

2 modes

6 customizable settings

2 handles

2 brush heads



HX8932/75



Removes up to 100% more stains in only 1 week*

Safely enhance your smile every day without giving up the staining foods and drinks you love! The Philips Sonicare whitening rechargeable toothbrush is proven to remove up to 100% more stains for whiter teeth in only 1 week*.

Proven to improve oral health

- Remove up to 7x more plaque with our InterCare brush head*
- Philips Sonicare technology helps brighten your smile
- Safe and gentle for your oral health needs

Optimize your brushing

- Choose from 2 modes and 3 intensity settings

Deep, gentle cleaning

- Sensitive Mode: Gentle teeth and gum cleaning

Provides a superior clean

- Clean and White Mode: proven to remove stains
- Angled neck to help tackle cleaning in tough areas
- Our unique technology brings you a powerful yet gentle clean

Designed around you

- Timers encourage thorough 2 minutes of brushing


Highlights

Gentle daily whitening



Naturally whiten teeth with our unique sonic technology. This brush head helps to remove everyday stains from food and drink so you can gently brighten your smile each time you brush.

Remove up to 7x more plaque



Click on our InterCare brush head to improve gum health in just 2 weeks. Extra-long bristles help remove more plaque in hard-to-reach spots and between teeth for healthy gums.

Clean and White Mode




2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*

Gentle teeth and gum cleaning



Philips Sonicare toothbrush with Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums

2 modes, 3 intensity settings



HealthyWhite+ lets you brush your way. Boost your brushing with 3 intensity settings, plus 2 modes to cover your whitening and cleaning needs. Choose Clean – for outstanding plaque removal in 2 minutes, and White – to remove stains like coffee or tea, and whiten teeth by up to 2 shades in only 2 weeks.

Philips Sonicare technology




Powerful sonic vibrations whip up your toothpaste into plaque fighting bubbles, and drive them deep between your teeth and along your gum line. At the same time your teeth will experience 62,000 gentle but effective brush movements. You'll get a month's worth of manual brushing in just 2 minutes.

Reach those tricky spots



The unique shape of your HealthyWhite+ toothbrush handle combined with our brush heads means that even tricky spots like back teeth get a thorough clean.

QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

Brush with confidence



Whatever your specific needs, you can be sure of a safe brushing experience with HealthyWhite+. Our sonic technology is suitable for use with braces, fillings, crowns and veneers, and can even be used in the daily treatment of gum disease.

Specifications

Modes

Clean: For exceptional everyday clean

White: Removes surface stains

3 intensities: High, Low, Medium

Items included

Handle: 2 HealthyWhite+

Brush heads: 2 InterCare standard

Travel case: 2

Charger: 2

Design and finishing

Color: Iridescent

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Removes up to 100% more plaque

Health benefits: Improves gum health in only two weeks

Whitening benefits: Whitens up to 2 shades in 2wks

Timer: Quadpacer and SmartTimer

Ease of use

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 3 weeks**

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

