

PHILIPS

sonicare

Sonic electric
toothbrush

DiamondClean

5 modes

2 brush heads

2 tongue brushes

Glass charger



HX9368/35



Whiter, healthier teeth for life

Best Philips Sonicare whitening sonic toothbrush

Sonicare's best whitening in our most elegant Philips Sonicare electric toothbrush. Switch to Sonicare.

Proven to improve oral health

- Up to 7x more plaque removal than a manual toothbrush
- A whiter smile in 1 week with our DiamondClean brush head*
- Improves gum health in just 2 weeks*
- Choose from 5 modes including Polish and Sensitive
- Safe and gentle for your oral health needs

Innovative technology

- A Philips Sonicare clean for your tongue

Provides a superior clean

- Philips Sonicare's advanced sonic technology
- Angled neck to help tackle cleaning in tough areas

Designed around you

- Timers encourage thorough 2 minutes of brushing
- Easy start program builds your Philips Sonicare routine



asimpleswitch.com


Highlights

Say goodbye to plaque




Densely-packed, high-quality bristles remove up to 7x more plaque than a manual toothbrush.

Up to 2x whiter teeth



Click on our DiamondClean brush head to gently yet effectively remove surface stains. Densely packed central stain removal bristles work hard to make your smile 2x whiter in just 7 days.*

Look after those gums



With optimal cleaning from your DiamondClean, gums will become healthier in 2 weeks*. Removing up to 7x more plaque along the gum line than a manual toothbrush*, you'll get your healthiest smile.

5 modes to clean your way



With DiamondClean, you get a refreshing clean every day. Our 5 modes cover all your brushing needs; Clean mode – for exceptional daily cleaning, Gum Care – to gently massage gums, Polish – to brighten your smile, Sensitive – for gentle yet effective cleaning of sensitive gums, and White – the ideal mode to remove surface stains.

Philips Sonicare technology




Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Sonic technology



TongueCare+ is the only tongue cleaner with Philips Sonicare strength, to clean away bad-breath bacteria with every stroke. 31,000 vibrations per minute break up tongue debris, sweep away bad-breath bacteria, and drive the bacteria-killing ingredients of the Philips Sonicare BreathRx tongue spray deeper. Your tongue will feel cleaner and your breath will stay fresh longer.

QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

Reach those tricky spots



The unique shape of your DiamondClean toothbrush handle combined with our brush heads means that even tricky spots like back teeth get a thorough clean.



[asimpleswitch.com](https://www.asimpleswitch.com)

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Modes

Clean: For exceptional everyday clean

Gum Care: Gently massages gums

Polish: Brighten and polish your teeth

Sensitive: Gentle teeth and gum cleaning

White: Removes surface stains

Items included

Handle: 2 DiamondClean

Brush heads: 2 DiamondClean standard

Glass charger: 1

Tongue brush: 2

Design and finishing

Color: Black, Pink

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Removes up to 7x more plaque*

Health benefits: Improves gum health in only two weeks

Whitening benefits: Whitens teeth 2X better*

Timer: Quadpacer and SmarTimer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Illuminated icon indicates battery life

Handle: Slim ergonomic design

Display: Illuminated display

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 3 weeks**

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

