

PHILIPS

sonicare

Sonic electric
toothbrush

For Kids

2 modes

1 brush head

10 stickers



HX6381/07



Sonic power maximizes your child's routine Electric toothbrush for kids

Philips Sonicare For Kids is an electric power toothbrush for kids aged 7 and older. It offers maximum plaque removal, sonic technology, customizable stickers and educational tools to help make proper brushing fun for a lifetime.

Designed to help build proper brushing habits

- KidTimer helps to increase brushing time
- 2 Kid-friendly power modes deliver a gentle, effective clean
- 2 brush head sizes available
- KidPacer encourage kids to brush effectively
- Fun interchangeable stickers for customization
- Multi grip design for parents and kids
- Preferred by dental professionals for their own children
- Anti-roll shape

Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Removes more plaque than a children's manual toothbrush

Safe for growing smiles

- Rubberized brush head is designed to protect young teeth

Highlights

Age-appropriate brush heads



This Philips Sonicare electric toothbrush has two brush head sizes that are specifically designed to gently clean and protect teeth as they grow

Power modes




With two kid-friendly power modes, this power toothbrush provides proper cleaning for different ages; a low mode for younger children and high mode for older kids.

Sonic technology



Features a unique dynamic action that gently and effectively reaches deep between teeth and along the gumline.

KidPacer




KidPacer encourage kids to brush their entire mouth effectively for the dentist-recommended 2 minutes

Stickers



Allows kids to make each brush uniquely their own

Preferred by dental profession




91% of dental professional parents prefer Philips Sonicare For Kids for their own children**

KidTimer



To help kids ease into a proper routine, this electric toothbrush slowly increases the brushing time over 90 days until reaching the dentist-recommended 2 minutes to help instill healthy habits naturally.

Maximum plaque removal



Philips Sonicare For Kids removes up to 75% more plaque than manual toothbrushes*

Rubberized brush head



The rubberized brush head of this electric toothbrush is designed to help protect growing teeth

Multi grip design



Ergonomically designed to help kids brush properly on their own

Specifications

Modes

Power modes: 2

Items included

Handles: 1 Sonicare for Kids

Brush heads: 1 Sonicare for Kids standard

Charger: 1

Stickers: 10 customization stickers

Design and finishing

Color: Aqua

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: 75% more effective*

Health benefits: For healthy oral care habits

Timer: KidTimer and Quadpacer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Light shows battery status

Handle: Slim ergonomic design, Rubber grip for easy handling

Display: Illuminated display

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 3 weeks**

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

