

PHILIPS

sonicare

Rechargeable sonic
toothbrush

EasyClean

1 mode

1 brush head

HX6511\02



Better plaque removal

Sonicare's unique dynamic cleaning action gently and effectively reaches deep between teeth and along the gumline.

Proven to improve oral health

- Clinically proven safe and gentle
- Helps whiten teeth

Offers you a customized cleaning experience

- Easy-start program to acclimate to the Sonicare experience

Provides a superior clean

- Dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- Removes up to 2x more plaque than a manual toothbrush
- Dynamic cleaning action drives fluid between teeth

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Highlights

Patented sonic technology



Sonicare's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck



Uniquely angled brush head neck makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Safe and gentle



Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

Smartimer




2-minute timer helps ensure dental professional recommended brushing time

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Specifications

Cleaning performance

Performance: Better plaque removal than a manual toothbrush

Health benefits: Helps improve gum health

Whitening benefits: Helps remove tea, coffee and other stains

Brush heads: 1 ProResults standard brush head for an All-Around Clean

Speed: Up to 31,000 brush strokes per minute

Timer: Smartimer and quadpacer

Ease of use

Battery indicator: Light in power button blinks to indicate when to recharge

Brush head system: Easy snap-on brush heads for optimal hygiene

Brushing time: 42 2-minute brushings OR 3 weeks

Handle: Slim ergonomic design for easy handling.

Items included

Brush head color rings: 3 color rings for easy sharing

Charger: 1 charger

Travel charger: 1

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Fully charged should provide 20 2-min. brushings

Design and finishing

Color: White

Service

Warranty: 2-year limited warranty

Power

Voltage: Multi-voltage charger

Features

Quadpacer: Encourages attentive brushing of each quadrant of the mouth

Rechargeable battery

Smartimer: 2-minute timer helps ensure recommended brushing

