

PHILIPS

sonicare

Sonic electric
toothbrush

HealthyWhite

2 modes

3 brush heads

HX6772/20



Naturally whitens teeth

The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Philips Sonicare HealthyWhite electric toothbrush is proven to remove up to 90% of stains for a whiter, brighter smile.

Proven to improve oral health

- Naturally whiter teeth with patented sonic technology
- Clinically proven safe and gentle

Provides a superior clean

- Removes up to 6 times more plaque than a manual toothbrush
- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

Offers a customized brushing experience

- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2 minute mode
- Easy-start program to acclimate to the Sonicare experience

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Highlights

Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

Clean Mode

For exceptional everyday clean. Helps maintain whitening treatments.

Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*

Removes up to 6x more plaque




Clinically proven to remove up to 6 times more plaque than a manual toothbrush.

Sonic technology



Our sonic technology creates a unique dynamic cleaning action that drives fluid deep between teeth and along the gumline for a thorough-yet-gentle clean

Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Smarter



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

Specifications

Modes

Clean: For exceptional everyday clean

Clean & White: 2 min Clean + 30 sec White

Items included

Handles: 1 HealthyWhite

Brush heads: 1 DiamondClean standard, 2 plaque control brush heads

Charger: 1

Travel case: 1

Design and finishing

Color: White

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Up to 100% more stain removal*

Health benefits: Helps improve gum health

Whitening benefits: Whitens teeth up to 2 shades

Timer: Quadpacer and SmarTimer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Light shows battery status

Handle: Slim ergonomic design, Rubber grip for easy handling

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 3 weeks**

Battery type: Lithium ION

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

