

PHILIPS

sonicare

Sonic electric
toothbrush

DailyClean 3100

1 mode

1 brush head

HX6201/57



Removes up to 7X more plaque*

*than a manual toothbrush

Offers you superior plaque removal at the push of a button

Proven to improve oral health

- Safe and gentle for your oral health needs

Innovative technology

- Philips Sonicare's advanced sonic technology

Designed around you

- Easy start program builds your Philips Sonicare routine
- Timers encourage thorough 2 minutes of brushing

Provides a superior clean

- Better plaque removal to help reduce cavities
- Angled neck to help tackle cleaning in tough areas

Highlights

Helps reduce cavities



Provides better plaque removal and access to back teeth; brushing twice a day every day with this brush can help reduce cavities

The only true sonic technology




Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Safe and gentle



Whatever your specific needs, you can be sure of a safe brushing experience with the 2 Series plaque control. Our sonic technology is suitable for use with braces, fillings, crowns and veneers, and can even be used in the daily treatment of gum disease.

Reach those tricky spots




The unique shape of your 2 Series toothbrush handle combined with our brush heads means that even tricky spots like back teeth get a thorough clean.

QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

Specifications

Modes

Clean: For exceptional everyday clean

Items included

Handle: 1 DailyClean 3100

Brush heads: 1 C2 Optimal Plaque Control

Charger: 1

Design and finishing

Color: Blue Roan

Cleaning performance

Performance: Removes up to 7x more plaque*

Health benefits: Helps reduce cavities

Timer: Quadpacer and Smartimer

Ease of use

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Brushing time: Up to 2 weeks

Technical specifications

Battery: Rechargeable

Battery type: NiMH

Operating time (full to empty): Up to 2 weeks

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

