



Sonic electric toothbrush

- 2 modes
- 1 brush head
- 8 Stickers

HX6311/07

Sonic power maximizes your child's routine

Electric toothbrush for kids

Philips Sonicare For Kids is an electric power toothbrush for kids aged 7 and older. It offers maximum plaque removal, sonic technology, customizable stickers and educational tools to help make proper brushing fun for a lifetime.

Benefits

Designed to help build proper brushing habits

- 2 brush head sizes available
- KidTimer helps to increase brushing time
- 2 Kid-friendly power modes deliver a gentle, effective clean
- Fun interchangeable K3 stickers for customization
- Multi grip design for parents and kids
- Preferred by dental professionals for their own children
- Anti-roll shape

Provides a superior clean

- Removes more plaque than a children's manual toothbrush
- Sonicare dynamic cleaning action drives fluid between teeth

Safe for growing smiles

- Rubberized brush head is designed to protect young teeth

PHILIPS

For Kids

Features

Age-appropriate brush heads



This Philips Sonicare electric toothbrush has two brush head sizes that are specifically designed to gently clean and protect teeth as they grow

KidTimer



To help kids ease into a proper routine, this electric toothbrush slowly increases the brushing time over 90 days until reaching the dentist-recommended 2 minutes to help instill healthy habits naturally.

Power modes



With two kid-friendly power modes, this power toothbrush provides proper cleaning for different ages; a low mode for younger children and high mode for older kids.

Stickers



Allows kids to make each brush uniquely their own.

Maximum plaque removal



Philips Sonicare For Kids removes up to 75% more plaque than manual toothbrushes*

Multi grip design



Ergonomically designed to help kids brush properly on their own

Sonic technology



Features a unique dynamic action that gently and effectively reaches deep between teeth and along the gumline.

Preferred by dental profession



91% of dental professional parents prefer Philips Sonicare For Kids for their own children**

Rubberized brush head



The rubberized brush head of this electric toothbrush is designed to help protect growing teeth

Anti-roll shape



Sturdy handle design allows kids to store upright as well as apply toothpaste while the toothbrush is laying flat.

Specifications

[field not found 'disclaimer_text']

Modes

Power modes 2

For Kids

Items included

Handles	1 Sonicare for Kids
Brush heads	1 Sonicare for Kids standard
Charger	1
Stickers	8 customization stickers 2 Bonus stickers

Design and finishing

Color	Aqua
-------	------

Cleaning performance

Speed	Up to 62000 brush movement/min
Performance	75% more effective*
Health benefits	For healthy oral care habits
Timer	KidTimer and Quadpacer

Ease of use

Brush head system	Easy click-on brush heads
Brushing time	Up to 3 weeks**
Battery indicator	Light shows battery status
Handle	Slim ergonomic design Rubber grip for easy handling
Display	Illuminated display

Technical specifications

Battery	Rechargeable
Battery type	Lithium ION
Operating time (full to empty)	Up to 3 weeks**

Power

Voltage	100-240 V
---------	-----------

Service

Warranty	2-year limited warranty
----------	-------------------------

Packaging dimensions

Height	24.40 cm
Width	11.10 cm
Depth	6.40 cm
Nett weight	0.30 kg
Gross weight	0.37 kg
EAN	0871010366878
Number of products included	1
Country of origin	CN
Harmonised Systems Code	850980

Outer Carton

Length	23.60 cm
Width	11.90 cm
Height	20.30 cm
Gross weight	1.25 kg
EAN	18710103668784
Number of consumer packagings	3

* Removes up to 7x more plaque vs. a manual toothbrush

* in hard to reach areas

** Philips Sonicare in home survey of U.S. dental professionals with children ages 4-10

*** based on two periods of two-minute brushings per day



data subject to change
2021, October 2

Version: 16.3.1
EAN: 08710103668787

© 2021 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com