

PHILIPS

sonicare

Sonic electric
toothbrush

3100 Series

Sonicare technology

Pressure sensor

QuadPacer and SmarTimer

Slim ergonomic design

HX3681/04



Goodbye manual toothbrush. Hello Sonicare.

Up to 3x better plaque removal*

Sonicare technology combined with our brush action gently removes plaque up to 3x better* than a manual toothbrush.

Designed around you

- 14-day long lasting battery life
- Ergonomic design makes the toothbrush easy to hold and use
- 2 intensity settings for a more powerful or gentle clean

Optimize your brushing

- Optimize your brushing with SmarTimer and QuadPacer

Easy transition to Sonicare

- Easy-Start to ease the transition

Protects your gums

- Built-in pressure sensor protects your teeth & gums

Provides superior clean

- Up to 3x better plaque removal vs a manual toothbrush*
- Our unique technology brings you a powerful yet gentle clean

3100 Series

Sonic electric toothbrush

HX3681/04

Highlights

Removes plaque 3x better*



The Sonicare electric toothbrush with advanced sonicare technology is clinically proven to remove plaque up to 3X better* than a manual toothbrush. It removes plaque from your teeth and along the gum line while protecting your gums.

Easy-Start




Our Easy-Start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

2 intensity settings



Switch between 2 intensity settings to make the brushing session comfortable for you.

Protects your gums



Built-in pressure sensor automatically detects pressure you apply, warns you, and reduces the toothbrush vibrations automatically to help protect your gums. The toothbrush will make a pulsing sound as a reminder to ease off the pressure. 7 out of 10 people found this feature helped them become better brushers.

SmartTimer and QuadPacer



2 minute SmartTimer and 30-second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

Advanced Sonicare Technology



Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two month's worth of manual brushing in just 2 minutes.** 31000 brush strokes per minute gently clean your teeth, break up plaque, and sweep it away for an exceptional daily clean.

Designed around you



Slim, light-weight, and ergonomic design makes the toothbrush easy to hold and use so that you can get an amazing clean without any strain.

Long lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

3100 Series

Sonic electric toothbrush

HX3681/04

Specifications

Power

Voltage: DC5V

Technical specifications

Operating time (full to empty): Up to 14 days

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: Black

Service

Warranty: 2-year limited warranty

Ease of use

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design

Items included

Charger: 1

Handle: 1 3000 Series

Brush head: 1 C1 ProResults standard

Cleaning performance

Performance: Removes up to 3x more plaque*

Timer: Quadpacer and SmarTimer

Pressure feedback: Vibrates handle to alert user

Modes

Clean: For exceptional everyday clean

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

