

PHILIPS

sonicare



Oral Irrigator

Power Flosser 3000

Quad Stream technology

Whisper quiet design

Pulse waves guide as you go



HX3711/20

Up to 180% more effective for healthier gums*

An effortless, thorough clean

Enjoy an effortless, thorough clean between your teeth and improve your gum health. Our Quad Stream technology helps you cover more area with less effort, while Pulse Wave technology gently guides you from tooth to tooth. Easy.

An effortless clean with advanced technology

- Easy to fill 550ml reservoir
- Quadstream technology: removes plaque from below the gumline
- Quad Stream technology for faster, more effective flossing
- Standard nozzle for a regular clean
- Pulse Wave technology guides from tooth to tooth
- Whisper quiet for great results without disturbing others
- Quad Stream nozzle for maximum coverage

Customize easily to your needs

- 2 flossing modes, 10 intensities


Power Flosser 3000

Oral Irrigator

HX3711/20


Highlights

Covers more area, effortlessly




X-shaped water streams cover more surface area for a faster, deeper clean versus manual floss. Effortlessly.

Cleans from below the gumline



The Quadstream technology removes up to 99% of plaque from 6mm deep periodontal pockets.

Guides from tooth to tooth



Gentle pulses of water guide you from tooth-to-tooth so that you don't miss a spot.

Whisper quiet design



Engineered for a discreet clean, so you get great results without disturbing others.

2 modes, 10 intensities




Enjoy a clean that's tailored to your needs. Clean mode uses a continuous flow of water for a great everyday clean, while Deep Clean uses distinct pulses for a more thorough clean. Whenever you need, adjust the intensity to your ideal comfort level.

F1 Standard nozzle



Use the standard nozzle for a regular clean. It uses a single stream of water to remove debris and clean thoroughly between your teeth.

F3 Quad Stream nozzle



For maximum coverage and an extra thorough clean, use the Quad Stream nozzle to concentrate water into four wide streams between teeth and along the gum line. The soft rubber guide is gentle on gums and every nozzle simply clicks on and off.

550ml reservoir



The 550ml reservoir holds enough water for one cleaning. To refill simply lift out, fill and click back in.

Power Flosser 3000

Oral Irrigator

HX3711/20

Specifications

Power

Voltage: Multi-voltage charger

Technical specifications

Energy Consumption: Off mode <0.3 W, Standby without display <0.5 W

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Nozzle attachment: Easily clicks on and off

Performance

Cleaning: Clean entire mouth in 60-90sec

Items included

F1 Standard nozzle: 1

F3 Quad Stream nozzle: 1

Power Flosser: 1

Modes

Clean: For exceptional everyday clean

Deep Clean+: For an invigorating deep clean

Intensities: 10

Health benefits

Gum health: Up to 180% more effective*

Plaque removal: Up to 3X as effective**

Sub-gingival: Removes up to 99% of plaque from below the gumline***



* vs. floss. When used with Quad Stream nozzle on setting 8 with a manual toothbrush in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss

** as floss in removing plaque between teeth and along the gumline when used with Quad Stream tip on setting 8 with a manual toothbrush in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss.

*** In a lab study (actual results may vary), in the deepest 2mm of a 6mm periodontal pocket model.