

PHILIPS

sonicare



Sonic electric toothbrush
with accessories - Pink

DiamondClean 9000

Connected brushing made easy

Built-in pressure sensor

Smart brush head recognition

4 modes, 3 intensities



HX9911/53

Whiter, healthier teeth for life

Whiter teeth in just 1 day*

Philips Sonicare's best whitening in our most elegant sonic electric toothbrush. Switch to Philips Sonicare.

Enjoy a stylish and user-friendly design

- Simple to charge and style to boot

DiamondClean users experience improved oral health

- Hard on plaque, gentle on your gums
- Remove up to 100% more stains* in just 3 days
- Removes up to 10 x more plaque* for a deep clean
- Stains disappear for a brighter, whiter smile

Let the DiamondClean 9000 adapt to your needs

- Personalised brushing experience
- Lets you know if you're pressing too hard
- BrushSync automatically selects the best mode for you
- Always know when to replace your brush heads

Stay on track with a personalised progress report

- Start and keep up healthy habits

DiamondClean 9000

Sonic electric toothbrush with accessories - Pink

HX9911/53


Highlights

Remove up to 5 x more stains***



Get rid of surface stains and enjoy a whiter smile with the W3 Premium White toothbrush head. The toothbrush head's densely packed bristles are especially designed to tackle stains and the brush is clinically proven to remove up to 100% more stains in just three days.***

Removes up to 10 x more plaque*



Know you're getting the deepest possible clean with the C3 Premium Plaque Control toothbrush head. The soft, flexible bristles are designed to curve around the contours of your teeth, giving you 4 x more surface contact** and up to 10 x more plaque removal from hard-to-reach spots.

W3 Premium White brush head




Click on the W3 Premium White toothbrush head to remove surface stains and reveal your whitest smile. Its densely packed central stain removal bristles whiten your teeth from day one.

C3 Premium Plaque Control



The C3 Premium Plaque Control toothbrush head is designed to give you your deepest clean yet. Soft, flexible sides and bristles perfectly curve around the contours of your teeth, giving you 4 x more surface contact and helping target hard-to-reach spots*.

Progress Report



The DiamondClean 9000 electric toothbrush gives you the guidance you need to improve and maintain healthy brushing habits in between

dentist check-ups. Built-in smart sensors let you know when you're using too much pressure, and by connecting your brushing experience to the Sonicare app, a personalised Progress Report helps you stay on track and see how much you've improved over time.

4 modes, 3 intensities



Philips Sonicare DiamondClean 9000 electric toothbrush comes with Clean, White+, Gum Health and Deep Clean+ to take care of your brushing needs. Clean mode is for exceptional everyday cleaning, White+ is for removing stains, Gum Health provides a gentle yet effective clean for gums and Deep Clean+ gives you an invigorating deep clean. Three intensities allow you to choose between a higher setting to boost your clean and a lower one for more sensitive mouths.

BrushSync mode pairing



Smart toothbrush heads ensure you're using the right mode and intensity for the best possible clean. For example, let's say you're using the W3 Premium White toothbrush head. Your DiamondClean 9000's BrushSync technology will automatically sync your toothbrush head with the White+ mode to help whiten your teeth.


DiamondClean 9000

Sonic electric toothbrush with accessories - Pink

HX9911/53

Highlights

Built-in pressure sensor



You might not notice if you're brushing too hard, but your Philips Sonicare DiamondClean 9000 will. If you need to ease up, the electric toothbrush will make a pulsing sound. It's a reminder to let your toothbrush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

BrushSync replacement reminder



All toothbrush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your toothbrush head for and how hard you've

been brushing. The BrushSync replacement reminder on your handle and a short beep remind you when it's time to replace it.

Long-life battery



Your Philips Sonicare DiamondClean 9000 electric toothbrush will last for up to two weeks on a single charge with regular use.

DiamondClean 9000

Sonic electric toothbrush with accessories - Pink

HX9911/53

Specifications

Product Update

Support Period: Philips will provide the necessary security updates until 31 December 2026.

Connectivity

Bluetooth® wireless technology: Connected brushing app

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): 14 days***

Battery: Rechargeable

Battery type: Lithium ION

Energy consumption: Standby without display <0.5 W

Design and finishing

Colour: Pink

Service

Warranty: 2-year limited warranty

Compatibility

Android compatibility: Android phones, Bluetooth 4.0-enabled tablets

iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS 7 operating system

Ease of use

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Replacement reminder: To always ensure best results, reminder icon lights up

Items included

Handle: 1 DiamondClean 9000

Brush heads: 1 C3 Premium Plaque Control

Travel case: USB charging travel case

Charger base: 1

Cleaning performance

Performance: Removes up to 10 x more plaque*

Timer: BrushPacer and SmarTimer

Whitening benefits: Whitens teeth in 1 day*

Pressure feedback: Vibrates handle to alert user

Modes

3 intensities: High, Low, Medium

Clean: For exceptional everyday cleaning

Deep Clean+: For an invigorating deep clean

Gum Health: Special attention to molars

White+: To remove surface stains

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

BrushSync Technology: smart brush head, Connects smart handle and

Software Support

Software updates: Philips offers relevant software updates for a period of 2 years after the date of purchase.

