

PHILIPS

AVENT

Manual Breast Pump & Store Set

Express, store and feed.

Pump with soft massage cushion

Includes 2 bottles and nipples

10x storage cups and lids

SCD221/10

 Product Picture

More comfort, more milk* Pump, store and feed

Manual breast pump set with storage cups

The Comfort Manual Breast Pump & Store set is ideal for moms who plan to express their breast milk. The convenient set includes all you need to pump and store your breast milk, feed your baby, and includes 2 bottles and 10 storage cups.

Easy to combine breast and bottle feeding

- Includes our Natural bottle and nipple for natural latch on

Easy to set up, express and clean

- Few separate parts and intuitive design
- Compact lightweight design

Clinically proven less leaning forward*

- Relaxing position, thanks to its special design

Store - Storage cups

- For safe storage and transport
- Pump in storage containers

Soft massage cushion with massaging petals

- Soft massage cushion gently stimulates milk flow

Highlights

Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, helps your milk to flow more easily.

Soft massage cushion



Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let down.

Easy on-the-go



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on-the-go more discreet.

Easy set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

Includes our Natural bottle



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Pump in storage containers



With the milk storage cup adapter, you can easily pump in the milk storage container. No leakage, no spillage of your precious breast milk.

Leak proof lid



Philips Avent storage cups have a lid that creates a secure seal for safe storage and transport.

Specifications

Design

Breast pump design: Compact design

Bottle design: Ergonomic shape, Wide neck

Material

Bottle: BPA free*, Polypropylene

Nipple: BPA free*, Silicone

Breast pump: BPA free* (food contact parts only)

Breast pads: Natural breathable material, Silky soft topsheet

Storage cups: BPA free*, Polypropylene

What is included

Breast pump body with handle: 1 pcs

Sealing disc for milk storage: 1 pcs

Milk storage cup adaptor: 2 pcs

Natural bottle 4oz: 1 pcs

Natural bottle 9oz: 1 pcs

Milk storage cup (180ml/6oz): 10 pcs

Breast pad sample packs: 1 pack (2 pcs)

Large massage cushion (25mm): 1 pcs (A smaller size cushion is available and sold separately)

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Development stages

Stage: 0 – 6 months

