

PHILIPS

sonicare



Tongue cleaning
starter kit

TongueCare+

1 tongue brush

1 tongue spray



HX8071/10

Clean tongue, fresh breath

Kills bad-breath bacteria

Click. Spray. Clean. Our breakthrough TongueCare+ tongue brush and spray instantly transforms your Philips Sonicare toothbrush into a powerful breath-freshening solution. Your tongue gets a Sonicare clean; you get long-lasting fresh breath.

An easy addition to your oral care routine

- From toothbrush to tongue cleaner in just one click

Our most advanced breath-freshening solution yet

- 240 MicroBristles for a gentle deep clean
- A Philips Sonicare clean for your tongue
- A triple threat to fight bad breath
- Kills bad-breath bacteria

Highlights

Triple clean



Up to 80% of bad breath comes from odour-producing bacteria that accumulate and hide deep within the porous surface of your tongue. TongueCare+ combines three powerful components that work together to form an advanced fresh-breath solution. The innovative tongue brush, with 240 MicroBristles, contours to the tongue to clean away bad breath bacteria below the surface. Specially formulated BreathRx tongue spray kills 99% of the germs that cause bad breath* and neutralises lingering odours. Sonicare sonic technology breaks up and cleans away debris with 31 000 strokes per minute.

Click-on tongue brush



With 240 rubber MicroBristles, the TongueCare+ tongue brush is specifically designed to clean the soft and porous surface of your tongue. The flexible MicroBristles contour to and gently clean around the ridges and grooves, to remove

bacteria build-up and drive our bacteria-killing BreathRx tongue spray deeper. The compact shape allows you to comfortably clean the entire tongue.

BreathRx tongue spray



The unique combination of ingredients in our Philips Sonicare BreathRx tongue spray is clinically proven to kill bad breath bacteria and neutralise odours instantly. The special thick formula coats the tongue and keeps it lubricated for a gentle cleaning experience. The cool mint flavour leaves your mouth feeling fresh and clean.


Easy to use



TongueCare+ tongue brushes quickly transform your Philips Sonicare toothbrush into a sonic-powered tongue cleaner. Each tongue brush clicks onto your Philips Sonicare toothbrush handle, just like a regular brush head. Simple to

use, replace and clean, these innovative tongue brushes are compatible with all Philips Sonicare click-on toothbrush handles, so it couldn't be easier to add tongue cleaning to your daily oral care routine. It fits all Philips Sonicare toothbrush handles except: PowerUp Battery and Essence.

Sonic power



TongueCare+ is the only tongue cleaner with Philips Sonicare strength, to clean away bad-breath bacteria with every stroke. 31 000 vibrations per minute break up tongue debris, sweep away bad-breath bacteria and drive the bacteria-killing ingredients of the BreathRx tongue spray deeper. Thanks to our Philips Sonicare technology, your tongue will feel cleaner and your breath will stay fresh for longer.

Specifications

Consumer Trade Item

Height: 22.5 cm

Width: 6.3 cm

Length: 4.1 cm

Net weight: 0.108 kg

Gross weight: 0.138 kg

GTIN: 8710103737698

Country of origin: US

Harmonized system code: 96032980

Outer Carton

Height: 26.51 cm

Width: 13.5 cm

Length: 13.5 cm

Net weight: 0.648 kg

Gross weight: 0.96 kg

GTIN: 18710103737695

Pallet

Height: 120.5 cm

Width: 100 cm

Length: 120 cm

Net weight: 145.152 kg

Gross weight: 247.322 kg

GTIN: 68710103737690

Pallet Layer

Height: 26.51 cm

Width: 100 cm

Length: 120 cm

Items included

Antibacterial tongue spray: 2 oz

Tongue brush: 1

Ingredients

Antibacterial tongue spray: 0.09% CPC, Zytex

Use by

Shelf life: 36 months

After opened: 6 months

