

The Philips logo, consisting of the word "PHILIPS" in a bold, blue, sans-serif font, is positioned within a white rectangular box that has a blue border on the right and bottom sides. This box is part of a larger blue header area.

# PHILIPS

Radiation Oncology



## Multimodality Simulation Workspace Platform

Philips Multimodality Simulation Workspace connects to Philips and non-Philips imaging devices or PACS to access image datasets including CT, MR, PET, Spectral CT, and Cone Beam CT. It provides vendor-neutral support for multimodality image fusion and contouring to help physicians better define tumor volume and surrounding organs-at-risk.

- Vendor-neutral, task-centered solution with easy, clear navigation
- Single space for multimodality image fusion and contouring
- Efficient, consistent simulation

# Minimum technical specifications

<b>Configurations available</b>	Single node
<b>Server</b>	Lenovo P520
<b>Hardware</b>	Intel(R) Xeon(R) W-2225 CPU @ 4.10GHz CPU(s): 8 Thread(s) per core: 2 Core(s) per socket: 4 RAM : 64 GB
<b>Operating systems</b>	Cent OS 7.9
<b>Network ports</b>	Qty 1 - 1Gbps
<b>Access client, MM Sim supported web browsers</b>	Google Chrome, and Edge v79 or later with JavaScript enabled to run OpenText Exceed TurboX v12
<b>Access client for MM Sim</b>	OpenText Exceed TurboX v12
<b>Access client OS requirements for MM Sim</b>	Windows 7 SP1, Windows 8.1 and Windows 10 64 bit, MacOS 10.14 (Mojave)   Display Resolution: 1280x1024 or greater
<b>Access client networking requirements</b>	10Mbps (minimum) 100Mbps or greater recommended.
<b>Remote serviceability</b>	Philips Remote Service Network
<b>Antivirus*</b>	McAfeeTP-10.7.1-45
<b>Security</b>	TLS 1.2 supported & Infra Hardening
<b>Security level (hardening)</b>	Cent OS (>80% of RMF benchmark)

\*Requires an active Service Maintenance Agreement

Multimodality Simulation Workspace is a new feature of the TumorLoc module, part of Pinnacle RTPS.

© 2021 Koninklijke Philips N.V. All rights are reserved. Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.



[www.philips.com](http://www.philips.com)  
4522 9917 2581 OCT 2021