




Philips Sonicare
ProtectiveClean 5100
Sonic electric toothbrush

Built-in pressure sensor

3 modes
1 x BrushSync feature
Travel case



HX6851/39

Whiter teeth. Gently does it.

Whitens teeth in just 1 week.

Feel the difference of a gentle clean with our Pressure Sensor while improving your gum health up to 100% more vs. a manual toothbrush

Proven to improve oral health

- Safe and gentle on sensitive areas, orthodontics and dental work
- Whitens teeth in just one week

Designed around you

- Encouragement to brush thoroughly
- Makes traveling easier for you

Innovative technology

- Lets you know if you're pressing too hard
- Always know when to replace your brush heads
- Philips Sonicare's advanced sonic technology
- Connects smart brush handle and smart brush heads

Optimise your brushing

- Choose from three modes

PHILIPS
sonicare

Sonic electric toothbrush
Built-in pressure sensor 3 modes, 1 x BrushSync feature, Travel case

HX6851/39


Highlights

Naturally whiter teeth




Click on the W2 Optimal White brush head to remove surface stains and reveal a whiter smile. With its densely packed central stain-removal bristles, it's clinically proven to whiten teeth in just one week.

Three modes



With a choice of three modes, this toothbrush allows you to customise brushing to your own needs. Clean mode is standard for superior cleaning. White is the ideal mode for removing surface stains. And Gum Care mode adds an extra minute of reduced-power brushing, so you can gently massage your gums.

Safe and gentle



You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns and veneers, and helps prevent cavities and improve gum health.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

Monitor your brushing pressure



Brushing too hard can damage your teeth and gums. To prevent this, your Philips Sonicare ProtectiveClean emits a gentle pulsing sound to remind you to ease off on the pressure.

BrushSync technology




A microchip-enabled technology that detects and synchronises the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

BrushSync replacement reminder



All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. When it's time to replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

Helpful timers



Need an electric toothbrush with a timer? Our QuadPacer lets you know when you've spent just the right amount of time cleaning each part of your mouth, while our Smartimer tells you when you've brushed for the recommended two minutes.

Handy travel case and charger



Our premium travel case lets you store your toothbrush hygienically, while our compact charging base keeps you topped up when you're on the go. You do get two weeks of regular use from a single full charge, but the charger's just the ticket for longer trips.

Sonic electric toothbrush
Built-in pressure sensor 3 modes, 1 x BrushSync feature, Travel case

HX6851/39

Specifications

Design and finishing

- Colour: Navy Blue

Power

- Voltage: 110-220 V

Technical specifications

- Operating time (full to empty): Up to 2 weeks
- Battery: Rechargeable
- Battery type: Lithium ION

Service

- Warranty: 2-year limited warranty

Ease of use

- Handle compatibility: Easy click-on brush heads
- Handle: Slim ergonomic design
- Battery indicator: Light shows battery status
- Brushing time: Up to 2 weeks

Items included

- Brush heads: 2 W2 Optimal White standard
- Travel case: 1
- Charger: 1
- Handle: 1 ProtectiveClean

Cleaning performance

- Health benefits: Improves gum health in only two weeks, Helps reduce cavities
- Speed: Up to 62,000 brush movements/min
- Timer: Quadpacer and Smartimer
- Pressure feedback: Vibrates handle to alert user

Modes

- Clean: For exceptional everyday cleaning
- Gum Care: Gently massages gums
- White: Removes surface stains

Smart sensor technology

- Pressure sensor: Alerts when brushing too hard



Issue date 2022-06-04

Version: 3.0.2

EAN: 08 71010 39126 20

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Removes up to 7x more plaque vs. a manual toothbrush