

PHILIPS

sonicare

Sonic electric
toothbrush

DailyClean 3100

1 mode

1 brush head

HX6221/56



Removes up to 2X more plaque*

*than a manual toothbrush

Offers you superior plaque removal at the push of a button

Designed around you

- Easy-start programme builds your Philips Sonicare routine
- Timers encourage thorough 2 minutes of brushing

Provides a superior clean

- Removes up to 2x more plaque than a manual toothbrush
- Better plaque removal to help reduce cavities

Optimise your brushing

- Philips Sonicare's advanced sonic technology
- Angled brush head for better access to back teeth

Highlights

Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

Helps reduce cavities


Provides better plaque removal and access to back teeth; brushing twice a day every day with the Philips Sonicare 3100 electric toothbrush can help reduce cavities

The only true sonic technology




Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

Better access to back teeth




This plaque-removing toothbrush has a slim, angled brush head which provides better access to back teeth and hard-to-reach places.

Quadpacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our Quadpacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

Ease into your new routine



Starting something new can take some getting used to. Our easy-start programme gives you the option of a gradual, gentle increase in brushing power with the first 14 uses of your new toothbrush.

Specifications

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: NiMH

Design and finishing

Colour: Blue Roan

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 2 weeks

Items included

Handle: 1 DailyClean 3100

Brush heads: 1 C1 ProResults standard

Charger: 1

Cleaning performance

Performance: Up to 2 X more plaque removal*

Speed: Up to 31,000 brush movements/min

Health benefits: Helps reduce cavities

Timer: Quadpacer and Smartimer

Modes

Clean: For exceptional everyday cleaning

