

# PHILIPS

## sonicare



Sonic electric  
toothbrush with app

ExpertClean 7400

Connected brushing made easy

Built-in pressure sensor

Smart brush head recognition

3 modes, 3 intensities



HX9685/03

## Improving patient routines for exceptional check-u

### Expert technology. Personalized guidance.

The Philips Sonicare ExpertClean focuses on how patients brush, helping them improve their oral care routine in between visits.

#### Proven to improve oral health

- Increased surface contact for eliminating plaque
- Better gum care for overall health

#### Improves patients' oral health

- Removes 10x more plaque\* for a deep clean
- Up to 7x healthier gums\* in just 2 weeks

#### Designed around your patient

- 3 modes, 3 intensities
- Alerts patients of low battery
- For the most effective cleaning, every day
- Guides patients to use correct amount of pressure
- Guides patients to better oral care habits
- Brush heads automatically select optimal settings


# ExpertClean 7400

Sonic electric toothbrush with app

HX9685/03

## Highlights

### Removes up to 10x more plaque\*




The C3 Premium Plaque Control brush head delivers our deepest clean. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact\*\* and removes up to 10x more plaque from hard-to-reach areas\*.

### Up to 7x healthier gums\*



The G3 Premium Gum Care brush head improves patient gum health. Smaller size and targeted gum line bristles deliver gentle, yet effective, cleaning along the gumline, where gum disease starts. It's clinically proven to deliver up to 100% less gum inflammation\* and up to 7x healthier gums in just two weeks.\*

### Progress Report



Sensors track and measure patients' brushing behaviors while the Philips Sonicare ExpertClean brush handle provides instant feedback. Over time, the Sonicare app retains and analyzes these habits and generates a Progress Report to help promote better brushing technique.

### 3 modes, 3 intensities



With three brushing modes (Clean, Gum Health, DeepClean+) and three intensity settings (Low, Medium, High), your patients can personalize their brushing experience for optimal comfort and a complete clean. See Specifications for more detailed information on modes.

### Smart Brush Head Recognition



A microchip in the brush head syncs with the Philips Sonicare ExpertClean handle. If a patient clicks on a C3 Premium Plaque

Control/Defense/Defence, G3 Premium Gum Care or W3 Premium White brush head, the toothbrush automatically selects the optimal brushing mode and intensity level. All your patients need to do is start brushing. It helps patients with having the best settings for the best performance and their needs

### Pressure sensor



If patients apply too much pressure, the intuitive pressure sensor will gently pulsate the handle and give patients the real-time feedback they need to adopt a softer touch. In a study of U.S. patients, 7 out of 10 aggressive brushers tested reduced their brushing pressure when using the pressure sensor.\*\*\*\*

### C3 Premium Plaque Control



C3 Premium Plaque Control is our deepest cleaning brush head yet. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact and remove up to 10x more plaque from hard-to-reach areas.\*

# ExpertClean 7400

Sonic electric toothbrush with app

HX9685/03

## Highlights

### G3 Premium Gum Care brush head



G3 Premium Gum Care is our best brush head for improving gum health. Its smaller size and targeted gumline bristles deliver gentle yet effective cleaning along the gumline, where gum disease starts. It's clinically proven to deliver up to 7x healthier gums\* in just 2 weeks.

### Brush head replacement



Brush heads exhibit fatigue and bristle wear after three months of normal use. Philips Sonicare ExpertClean provides an accurate smart brush head replacement reminder based on how often patients brush, how long they brush and how much pressure they use. When it's time for a new brush head, a light on the handle lets patients know, so they can continue to achieve their best possible results.

### Battery indicator



When fully charged, Philips Sonicare ExpertClean has enough battery for three weeks of regular use. The battery indicator illuminates and clearly alerts patients when the battery is low and requires charging.

# ExpertClean 7400

Sonic electric toothbrush with app

HX9685/03

## Specifications

### Connectivity

**Bluetooth® wireless technology:** Connected brushing app

### Power

**Voltage:** 110-220 V

### Technical specifications

**Operating time (full to empty):** 14 days\*\*\*

**Battery:** Rechargeable

**Battery type:** Lithium ION

### Design and finishing

**Color:** White and Light Blue

### Service

**Warranty:** 2-year limited warranty

### Compatibility

**Android compatibility:** Android phones, Bluetooth 4.0 enabled tablets

**iOS compatibility:** iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

### Ease of use

**Handle:** Slim ergonomic design

**Handle compatibility:** Easy click-on brush heads

**Replacement reminder:** To always ensure best results, reminder icon lights up

### Items included

**Handle:** 1 ExpertClean

**Travel case:** 1

**Brush heads:** 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care, 1 W3 Premium White

**Charger:** 1

### Cleaning performance

**Health benefits:** Up to 7x healthier gums\*

**Performance:** Removes up to 10x more plaque\*

**Whitening benefits:** Helps naturally remove stains

**Timer:** BrushPacer and SmarTimer

**Pressure feedback:** Vibrates handle to alert user

### Modes

**Clean:** For exceptional everyday clean

**3 intensities:** Low, Medium, High

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** Special attention to molars

### BrushSync Mode Pairing

**C3 Premium Plaque Control:** Pairs with Clean mode

**G3 Premium Gum Care brush head:** Pairs with Gum Health mode

### Smart sensor technology

**BrushSync Replacement Reminder:** Always know when to, replace brush heads

**Pressure sensor:** Alerts when brushing too hard

**BrushSync Technology:** Connects smart handle and, smart brush head

**Location sensor:** Tracks & improves coverage

**Scrubbing sensor:** Guides to reduce scrubbing

