

# PHILIPS

## sonicare

Sonic electric  
toothbrush

FlexCare

3 modes

1 brush head

HX6911/76



## Advanced Cleaning

Toothbrush for gentle cleaning of teeth and gum

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare electric toothbrush that adapts to your oral care needs.

### **Proven to improve oral health**

- Improves gum health in only two weeks

### **Guides you to follow dental recommendations**

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

### **Offers a customized brushing experience**

- Easy-start program to acclimate to the Sonicare experience
- Three different modes for a better cleaning experience

### **Provides a better clean in between teeth**

- Safe on orthodontics, dental work and veneers
- Sonicare dynamic cleaning action drives fluid between teeth


# Highlights

## Improves gum health




This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

## Three different cleaning modes



The three modes include a Clean mode for an outstanding clean in 2 minutes. A White mode for removing surface stains and a Massage mode to stimulate the gum.

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## SmarterTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Safe and gentle



The gentle motion of this Philips Sonicare power toothbrush is safe for dental work, including orthodontics, implants and veneers, so you can brush without worry.

## Sonic technology



Our sonic technology creates a unique dynamic cleaning action that drives fluid deep between teeth and along the gumline for a thorough-yet-gentle clean.



# Specifications

## Power

**Voltage:** 110–220 V

## Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** Up to 2 weeks

**Battery type:** Lithium ION

## Design and finishing

**Color:** Black

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Battery indicator:** Light shows battery status

**Brush head system:** Easy click-on brush heads

**Handle:** Slim ergonomic design

**Brushing time:** Up to 2 weeks

## Items included

**Charger:** 1

**Handles:** 1 FlexCare

**Travel case:** 1

**Brush heads:** 1 ProResults standard

## Cleaning performance

**Health benefits:** Helps improve gum health

**Speed:** Up to 62000 brush movement/min

**Timer:** Quadpacer and SmarTimer

## Modes

**Clean:** For exceptional everyday clean

**White:** Removes surface stains

**Massage:** Invigorating massage

**2 Routines:** Go Care and Max Care

