

# PHILIPS

## sonicare

Sonic electric  
toothbrush - Dispense

EasyClean

1 mode

2 brush heads

HX6582/51



## Healthier gums in 2 weeks

The switch from a manual toothbrush made easy

Philips Sonicare EasyClean combines our patented Sonicare toothbrush technology and the ProResults gum health brush head to remove up to 6x more plaque\*. It offers an effective solution to prevent and reduce the early signs of gum disease.

### Designed around you

- EasyStart
- Recharge indicator

### Guides you to follow dental recommendations

- Quadpacer
- Smartimer

### Provides a superior clean

- Designed for optimal gum health improvement
- Patented Sonicare toothbrush technology
- Exclusive bonus items only available to dental professionals

## Highlights

### Bonus items not sold at retail



Now you can offer an exclusive Philips Sonicare EasyClean bonus item your patients won't find anywhere else, for improved gum care. Designed around patient needs for dental offices like yours, the package includes: 1 extra ProResults gum health brush head.

### A gentler experience




This Philips Sonicare toothbrush head features densely-packed, high-quality tufts to effectively remove plaque and improve gum health. Specially curved power tufts to easily reach teeth at the back of the mouth. Also features gentle trim profile to cushion teeth for an extra gentle brushing experience.

### Sonic technology




Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

### EasyStart



Gently increases power over the first 14 uses to ease patients into the Sonicare brushing experience.

### Quadpacer



Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

### Smartimer



Two-minute timer helps ensure dental professional recommended brushing time.

### Recharge indicator

Light blinks to let user know the battery needs recharging.

## Specifications

### Power

Voltage: 110-220 V

### Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: Lithium ION

### Design and finishing

Color: Glacier Green

### Service

Warranty: 2-year limited warranty

### Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 2 weeks

### Items included

Handles: 1 EasyClean

Brush heads: 2 ProResults gum health

Charger: 1

### Cleaning performance

Health benefits: Improve gum health in 2 weeks\*

Performance: Removes up to 6X more plaque\*

Whitening benefits: Helps maintain white teeth

Timer: Quadpacer and Smartimer

Speed: Up to 62000 brush movement/min

### Modes

Clean: For exceptional everyday clean

