

# PHILIPS

## sonicare

Sonic electric  
toothbrush

DailyClean 3100

1 mode

1 brush head

HX6221/67



## Removes up to 2X more plaque\*

\*than a manual toothbrush

Offers you superior plaque removal at the push of a button

### Designed around you

- Easy-start programme builds your Philips Sonicare routine
- Timers encourage thorough 2 minutes of brushing

### Provides a superior clean

- Removes up to 2x more plaque than a manual toothbrush
- Better plaque removal to help reduce cavities

### Optimise your brushing

- Philips Sonicare's advanced sonic technology
- Angled brush head for better access to back teeth

# Highlights

## Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

## Helps reduce cavities


Provides better plaque removal and access to back teeth; brushing twice a day every day with the Philips Sonicare 3100 electric toothbrush can help reduce cavities

## The only true sonic technology




Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

## Better access to back teeth




This plaque-removing toothbrush has a slim, angled brush head which provides better access to back teeth and hard-to-reach places.

## Quadpacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our Quadpacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

## Ease into your new routine



Starting something new can take some getting used to. Our easy-start programme gives you the option of a gradual, gentle increase in brushing power with the first 14 uses of your new toothbrush.

# Specifications

## Power

Voltage: 110-220 V

## Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: NiMH

## Consumer Trade Item

Height: 23.7 cm

Width: 17.3 cm

Length: 6.5 cm

Net weight: 0.27 kg

Gross weight: 0.401 kg

GTIN: 8710103912606

Country of origin: CN

Harmonized system code: 85098000

## Outer Carton

Height: 25.6 cm

Width: 18.65 cm

Length: 21.4 cm

Net weight: 0.81 kg

Gross weight: 1.402 kg

GTIN: 18710103912603

## Design and finishing

Colour: Black and White

## Service

Warranty: 2-year limited warranty

## Pallet Layer

Height: 25.6 cm

Width: 100 cm

Length: 120 cm

## Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 2 weeks

## Items included

Charger: 1

Brush heads: 1 C1 ProResults standard

Handle: 1 DailyClean 3100

## Cleaning performance

Performance: Up to 2 X more plaque removal\*

Speed: Up to 31,000 brush movements/min

Health benefits: Helps reduce cavities

Timer: Quadpacer and Smartimer

## Modes

Clean: For exceptional everyday cleaning

