

PHILIPS

sonicare

Sonic electric
toothbrush

ProtectiveClean
4700

Built-in pressure sensor

2 modes, 3 intensities

1 x BrushSync feature

HX6481/11



Cleaner teeth. Gently does it.

Remove up to 10x more plaque than a manual

Feel the difference of a gentle clean with our Pressure Sensor while removing up to 10x more plaque than a manual toothbrush

Proven to improve oral health

- Safe & gentle on sensitive areas, orthodontics & dental work

Designed around you

- Encouragement to brush thoroughly
- Always know when to replace your brush heads

Provides a superior clean

- Up to 10x more plaque removal than a manual toothbrush

Innovative technology


- Lets you know if you're pressing too hard
- Connects smart brush handle and smart brush heads
- Philips Sonicare's advanced sonic technology

Optimize your brushing

- Two modes, three intensity settings


Highlights

Say goodbye to plaque



And thanks to its flexible design, it removes up to 10x more plaque from hard-to-reach areas* for a truly deep clean along the gum line and between teeth.

Brushing your way



This toothbrush allows you to customize your brushing with a choice of two modes and three intensities. Clean mode is standard for superior cleaning. And Gum Care mode adds an extra minute of reduced-power brushing, so you can gently massage your gums. Three intensities allow you to switch the order of 'higher' and 'lower' to align with the intensities from left to right.

Safe and gentle



You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns, and veneers, and helps prevent cavities and improve gum health.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Monitor your brushing pressure



Brushing too hard can damage your teeth and gums. To prevent this, your Philips Sonicare ProtectiveClean emits a gentle pulsing sound to remind you to ease off on the pressure.

BrushSync replacement reminder




All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. When it's time to replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

BrushSync technology



A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

Helpful timers



Need an electric toothbrush with a timer? Our QuadPacer lets you know when you've spent just the right amount of time cleaning each part of your mouth, while our Smartimer tells you when you've brushed for the recommended two minutes.

Specifications

Power

Voltage: 110–220 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 2 weeks

Battery type: Lithium ION

Design and finishing

Color: White and Mint

Service

Warranty: 2-year limited warranty

Ease of use

Battery indicator: Light shows battery status

Handle compatibility: Easy click-on brush heads

Handle: Slim ergonomic design

Brushing time: Up to 2 weeks

Replacement reminder: To always ensure best results, reminder icon lights up

Items included

Charger: 1

Brush heads: 1 C3 Premium Plaque Control

Handle: 1 ProtectiveClean

Cleaning performance

Performance: Removes up to 10x more plaque*

Health benefits: Helps improve gum health,

Helps reduce cavities

Timer: Quadpacer and SmarTimer

Pressure feedback: Vibrates handle to alert user

Modes

Clean: For exceptional everyday clean

3 intensities: Low, Medium, High

Gum Care: Gently massages gums

Smart sensor technology

BrushSync Replacement Reminder: Always

know when to, replace brush heads

Pressure sensor: Alerts when brushing too hard

