

# PHILIPS

## sonicare

Sonic electric  
toothbrush

DailyClean 2100

1 mode

1 brush head



HX3211/33



Clinically proven to provide superior cleaning\*  
Improve gum health up to 100%\*

Designed to be the perfect toothbrush to transition from your manual toothbrush. Until now, it feels like you've never brushed your teeth before.

#### Proven to improve oral health

- Improves gum health by up to 100%\*

#### Designed around you

- Timers encourage thorough 2 minutes of brushing
- Easy start program builds your Philips Sonicare routine

#### Innovative technology

- Philips Sonicare's advanced sonic technology

#### Optimize your brushing

- Battery life
- Click-on brush head system

# Highlights

## The only true sonic technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

## Improves gum health



Densely-packed, high-quality bristles give you an extra gentle brushing to remove plaque along the gumline to improve gum health up to 100% more than a manual toothbrush. And, the specially curved power tip makes reaching the teeth at the back of your mouth a breeze.

## QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our


Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

## Battery life




Lasts up to 14 days between charges.

## Click-on brush head system



This brush head clicks on and off your brush handle for a secure fit and easy maintenance and cleaning. It fits Philips Sonicare toothbrush PowerUp Battery, Essence+/Elite+/CleanCare+ handles.

## Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

# Specifications

## Power

Voltage: 110-220 V

## Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 14 days

Battery type: NiMH

## Design and finishing

Color: Dark blue

## Service

Warranty: 2-year limited warranty

## Ease of use

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Handle: Ergonomic design, Rubber grip for easy handling

Brushing time: Up to 10 days

## Items included

Charger: 1

Brush heads: 1 G2 Optimal Gum Care

Handle: 1 DailyClean 2100

## Cleaning performance

Health benefits: Improves gum health

Timer: Quadpacer and Smartimer

## Modes

Clean: For exceptional everyday clean

