



Philips Sonicare
HealthyWhite
Sonic electric toothbrush

3 modes
2 brush heads

HX6732/74



Whiter, healthier teeth

The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Sonicare HealthyWhite electric toothbrush HX6711/02 is proven to remove everyday stains in just two weeks, with the use of the Clean & White Mode.

Proven to improve oral health

- Philips Sonicare toothbrush helps whiten teeth

Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

Provides a superior clean

- Improves gum health in just 2 weeks
- Sonicare dynamic cleaning action drives fluid between teeth
- Removes up to 7x more plaque than a manual toothbrush

Offers a customized brushing experience


- Easy-start program to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains

PHILIPS

sonicare

Highlights

Removes up to 7x more plaque



Removes up to 7 times more plaque from hard to reach places than a manual toothbrush.

Extra-long bristles




The InterCare brush head features extra-long, densely-packed, high-quality bristles to target hidden plaque caught deep between your teeth. It is clinically proven to reduce gingival bleeding and inflammation in just 2 weeks.

Helps whiten teeth




This Philips Sonicare electric toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Clean and White Mode




2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*

Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth

and signals for you to move on resulting in more consistent cleaning throughout your mouth

SmarterTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.



Sonic electric toothbrush
3 modes 2 brush heads

HX6732/74

Specifications

Design and finishing

- Color: White

Power

- Voltage: 110-220 V

Technical specifications

- Operating time (full to empty): Up to 3 weeks**
- Battery: Rechargeable
- Battery type: Lithium ION

Service

- Warranty: 2-year limited warranty

Ease of use

- Battery indicator: Light shows battery status
- Brush head system: Easy click-on brush heads
- Handle: Slim ergonomic design, Rubber grip for easy handling
- Brushing time: Up to 3 weeks**

Items included

- Charger: 2
- Handles: 2 HealthyWhite
- Travel case: 2
- Brush heads: 2 InterCare standard

Cleaning performance

- Performance: Removes up to 7x more plaque*
- Health benefits: Improves gum health in only two weeks
- Speed: Up to 62000 brush movement/min
- Timer: Quadpacer and SmarTimer
- Whitening benefits: Helps whiten teeth

Modes

- Clean: For exceptional everyday clean
- Sensitive: Gentle teeth and gum cleaning
- Clean & White: 2 min Clean + 30 sec White



Issue date 2021-12-27

Version: 3.3.2

EAN: 00 07502 00461 54

© 2021 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Removes up to 7x more plaque vs. a manual toothbrush
**based on two periods of two-minute brushings per day, on clean mode