



Philips Sonicare C3
Premium Plaque Control
Standard sonic toothbrush
heads

3-pack
Standard size
Click-on
BrushSync mode pairing

HX9043/25

Our deepest cleaning brush head

Up to 10x more plaque removal* for great results

Patients can achieve their deepest clean yet. With its soft flexing sides, our Premium plaque control brush head adapts to the contours of teeth and gums to provide 4x more surface contact than a DiamondClean brush head for a complete clean

Innovative technology

- Handle automatically selects optimal mode and intensity***
- Philips Sonicare's advanced sonic technology

Designed to optimize performance

- Rigorously tested to meet patients' oral health needs
- Works with any Philips Sonicare click-on toothbrush
- For the most effective cleaning, every day

Philips Sonicare's deep clean

- Up to 4x more surface contact** for an effortless deep clean

For exceptional plaque removal

- Up to 10x more plaque removal*, even in hard-to-reach areas

PHILIPS

sonicare


Highlights

Say goodbye to plaque




The Philips Sonicare C3 Premium Plaque Control brush head gives you our deepest clean. Thanks to its flexible sides and unique bristle pattern design, it removes up to 10x more plaque, even in hard-to-reach areas* for a truly deep clean along the gum line and between teeth.

Philips Sonicare technology



Patients get a personalized clean every time they brush with our adaptive cleaning technology. Soft, flexible rubber sides let C3 Premium Plaque Control adapt to the unique contours of the mouth. Our bristles adjust to gums and teeth, giving patients up to 4x more surface contact than a DiamondClean brush head for deeper cleaning, even in hard-to-reach spots. Adaptive cleaning technology also allows gentle tracking along the gumline; absorbs excessive brushing pressure; and enables an enhanced sweeping motion, for a unique in-mouth feel and superior cleaning.

BrushSync mode pairing




Patients will always get their best possible clean with our BrushSync mode pairing feature. A microchip in the C3 Premium Plaque Control brush head syncs with the connected Philips Sonicare toothbrush handle. It prompts the handle to automatically select the optimal brushing mode and intensity level for exceptional plaque removal. All your patients need to do is start brushing.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Smart brush head replacement



Brush heads exhibit fatigue and bristle wear after three months of normal use. Philips


Sonicare smart toothbrushes provide accurate brush head replacement reminders based on how often patients brush and how much pressure they use. For those who don't use a smart toothbrush, C3 Premium Plaque Control's blue reminder bristles fade to white to let patients know it's time for a replacement brush head.

Easy to click-on



Premium White simply clicks on to a Philips Sonicare toothbrush handle for a secure fit and easy maintenance and cleaning. Premium White fits all Philips Sonicare toothbrush handles except PowerUp Battery and Essence.

Safe on teeth and gums



All genuine Philips Sonicare brush heads are proven safe and gentle on teeth and gums. Our products are strictly tested to ensure patients get exceptional performance and durability each time they brush.

Specifications

Design and finishing

- Color: Black

Design and finishing

- Bristle stiffness feel: Medium soft
- Smart Brush Head Recognition
- Material brush head: Soft, flexible rubber sides
- Reminder bristles: Blue bristle color fade away
- Size: Standard

Compatibility

- Brush head system: Click-on
- Suitable for these models: 2 Series plaque control, 2 Series plaque defense, 3 Series gum health, DiamondClean, DiamondClean Smart, EasyClean,

Essence+, FlexCare, FlexCare Platinum, FlexCare Platinum Connected, FlexCare+, for Kids, HealthyWhite, HealthyWhite+, PowerUp, ProtectiveClean

Items included

- Brush heads: 3 C3 Premium Plaque Control

Quality and performance

- Replacement: Every 3 months
- Tested: for optimal usage

Health benefits

- Plaque removal: Removes 10x more plaque*



Issue date 2022-07-30

Version: 3.3.1

EAN: 00 07502 00648 75

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Removes up to 7x more plaque vs. a manual toothbrush
* *than a DiamondClean brush head
* **BrushSync™ mode-pairing is only compatible with Philips Sonicare
BrushSync™ enabled toothbrush handles